

**ABSTRACT**

**Background:** Cefovecin (FOV) (Pfizer Animal Health, Kalamazoo, MI), is a new semi-synthetic extended-spectrum cephalosporin antibiotic. It is indicated for treatment of bacterial infections in small animals. An evaluation was undertaken to determine the accuracy and reproducibility of the Sensititre 18 to 24-hour dried susceptibility system (manufactured by TREK Diagnostic Systems, Cleveland, OH) with Cefovecin, compared to the CLSI (M31) reference broth microdilution method (BMD). The range tested for Cefovecin was 0.0005 – 32 µg/ml on both the Sensititre and BMD plates.

**Materials and Methods:** Two hundred isolates (20 *Staphylococcus intermedius*, 20 *Staphylococcus aureus*, 20 *Streptococcus canis* (group G), 20 *Escherichia coli*, 20 *Proteus mirabilis*, 20 *Pasteurella multocida*, 20 *Histophilus somnus*, 20 *Actinobacillus pleuropneumoniae*, 20 *Bacteroides fragilis* and 20 *Prevotella* species) were tested with FOV and 3 comparators, Cefpodoxime, Ceftiofur, and Cefoxitin, using the Sensititre 18-24 h dried susceptibility plate and by CLSI (M31) BMD for comparison. For reproducibility testing, 10 isolates were tested 3 times daily over a period of three days. Recommended CLSI quality control (QC) organisms were tested daily and were within the CLSI expected quality control ranges for FOV and all comparators.

**Results:** Comparison of the Sensititre plate to the CLSI (M31) BMD resulted in 100% essential agreement (+/- one log<sub>2</sub> dilution) for FOV and 100% for comparators. Reproducibility was calculated as the percentage of results within +/- one log<sub>2</sub> dilution of the modal value. Overall agreement for the reproducibility of FOV and for comparators was 100%.

**Conclusions:** This evaluation indicates that the Sensititre 18 to 24-hour dried susceptibility system with FOV is equivalent to the CLSI (M31) BMD and is an acceptable method for susceptibility testing of FOV.

**PURPOSE OF THE STUDY**

To evaluate the performance of Cefovecin, an extended-spectrum cephalosporin antibiotic, on the Sensititre 18 – 24 hour susceptibility system compared to the CLSI reference broth microdilution method (M31).

**MATERIALS & METHODS**

**Organisms Tested**

- 200 Clinical and Challenge isolates consisting of:

Organism	Number of Isolates Tested
<i>Staphylococcus intermedius</i>	20
<i>Staphylococcus aureus</i>	20
<i>Streptococcus canis</i> (grp G)	20
<i>Escherichia coli</i>	20
<i>Proteus mirabilis</i>	20
<i>Pasteurella multocida</i>	20
<i>Histophilus somnus</i>	20
<i>Actinobacillus pleuropneumoniae</i>	20
<i>Bacteroides fragilis</i> (group)	20
<i>Prevotella</i> species	20

- 10 Reproducibility isolates with on-scale endpoints

Organisms Tested
<i>Actinobacillus pleuropneumoniae</i>
<i>Actinobacillus pleuropneumoniae</i> ATCC 27090
<i>Streptococcus canis</i>
<i>Streptococcus pneumoniae</i> ATCC 49619
<i>Streptococcus intermedius</i>
<i>Bacteroides fragilis</i> ATCC 25285
<i>Prevotella</i> species
<i>Staphylococcus aureus</i> ATCC 29213
<i>Escherichia coli</i>
<i>Escherichia coli</i> ATCC 25922

**MATERIALS & METHODS cont.**

- CLSI quality control strains (20 replicates of each tested)

Organism	ATCC Number
<i>Staphylococcus aureus</i>	ATCC 29213
<i>Histophilus somnus</i>	ATCC 700025
<i>Actinobacillus pleuropneumoniae</i>	ATCC 27090
<i>Streptococcus pneumoniae</i>	ATCC 49619
<i>Escherichia coli</i>	ATCC 25922
<i>Eubacterium lentum</i>	ATCC 43055
<i>Bacteroides fragilis</i>	ATCC 25285
<i>Bacteriodes thetaiotaomicron</i>	ATCC 29741

\*Daily Quality Control isolates were run on both the reference and dried plates each test day.

**Antimicrobials Tested**


Antimicrobials Tested	Range Tested µg/ml	Supplied By
Cefovecin	0.0005-32	Pfizer Animal Health Kalamazoo, MI
Cefpodoxime	0.015-32	Sankyo, Japan
Ceftiofur	0.00025-0.25	Pfizer
Cefoxitin	0.25-256	Groton, CT Sigma-Aldrich Saint Louis, MO

**SUSCEPTIBILITY TESTING METHODS**

- Each isolate was tested using a custom dried Sensititre 18 – 24h susceptibility plate containing Cefovecin and comparator drugs. The plates were set up according to manufacturers' instructions.
- The CLSI reference broth microdilution plates were prepared and tested according to the Clinical Laboratory Standards Institute (CLSI- M31).
- Repeat testing was performed if any organism/antimicrobial combination resulted in a greater than +/- one log<sub>2</sub> dilution error. These would be repeated in triplicate on both the dried test plate and the reference plate.
- Reproducibility: 10 Clinical/ challenge strains that had on-scale endpoints, (chosen from the 200 previously chosen isolates) including the CLSI ATCC QC isolates, were tested in triplicate on 3 separate days. Any results on reproducibility with a greater than +/- one log<sub>2</sub> dilution difference would also be re-tested as per the criteria for repeat testing.


**RESULTS**

***Staphylococcus aureus, Staphylococcus intermedius***




**RESULTS cont.**


***Histophilus somnus, Actinobacillus pleuropneumoniae***




***Bacteriodes fragilis, Prevotella species***



***Escherichia coli, Proteus mirabilis***



***Streptococcus canis, Pasteurella multocida***



**RESULTS cont.**

**Reproducibility**

Isolate	Day 1 TEST FOV			Day 2 TEST FOV			Day 3 TEST FOV			MODE
	A	B	C	A	B	C	A	B	C	
<i>A. pleuropneumoniae</i>	0.015	0.03	0.03	0.03	0.03	0.03	0.03	0.03	0.03	0.03
<i>A. pleuropneumoniae</i> 27090	0.015	0.015	0.015	0.015	0.03	0.03	0.015	0.03	0.015	0.015
<i>S. canis</i>	0.06	0.06	0.06	0.06	0.06	0.06	0.06	0.06	0.06	0.06
<i>S. pneumoniae</i> 49619	0.25	0.25	0.25	0.25	0.25	0.25	0.12	0.12	0.12	0.25
<i>S. intermedius</i>	1	1	1	1	1	1	1	1	1	1
<i>S. aureus</i> 29213	1	1	1	1	1	1	1	1	1	1
<i>E. coli</i>	1	1	1	1	1	1	1	1	2	1
<i>E. coli</i> 25922	1	1	1	1	1	1	1	1	1	1
<i>Prevotella</i> spp.	0.25	0.25	0.25	0.25	0.25	0.25	0.25	0.25	0.25	0.25
<i>B. fragilis</i> 25285	8	16	8	8	8	8	16	8	8	8

Isolate	Day 1 REFERENCE FOV			Day 2 REFERENCE FOV			Day 3 REFERENCE FOV			MODE
	A	B	C	A	B	C	A	B	C	
<i>A. pleuropneumoniae</i>	0.015	0.015	0.015	0.015	0.015	0.015	0.015	0.015	0.015	0.015
<i>A. pleuropneumoniae</i> 27090	0.015	0.015	0.015	0.015	0.015	0.015	0.015	0.015	0.015	0.015
<i>S. canis</i>	<-0.03	<-0.03	<-0.03	0.06	0.06	0.06	0.06	0.06	0.06	0.06
<i>S. pneumoniae</i> 49619	0.12	0.12	0.12	0.12	0.12	0.12	0.12	0.12	0.12	0.12
<i>S. intermedius</i>	1	1	1	1	1	1	0.5	1	1	1
<i>S. aureus</i> 29213	0.5	0.5	0.5	0.5	0.5	0.5	0.5	0.5	0.5	0.5
<i>E. coli</i>	1	1	1	1	1	1	1	1	1	1
<i>E. coli</i> 25922	0.5	1	0.5	0.5	0.5	0.5	0.5	0.5	0.5	0.5
<i>Prevotella</i> spp.	0.25	0.25	0.25	0.25	0.25	0.25	0.25	0.25	<-0.12	0.25
<i>B. fragilis</i> 25285	8	8	8	8	8	8	8	8	8	8

**CONCLUSION**

The Sensititre 18-24 hour dried susceptibility system demonstrated an equivalent level of performance to the CLSI, M31 reference broth microdilution plate for susceptibility testing of Cefovecin.

**Isolates Tested:**

- The overall essential agreement for Cefovecin, within a +/- one log<sub>2</sub> dilution range, was 100% using the manual method.
- The overall essential agreement for comparators, within a +/- one log<sub>2</sub> dilution range, was 100% using the manual method.

Essential agreement for the 200 isolates were calculated using the +/- one log<sub>2</sub> dilution standard for comparison studies. The agreement rates for the 200 isolates were as follows:

Antimicrobials	No. Isolates Tested	% Essential Agreement
Cefovecin		