

## Multi-Site Evaluation of Moxifloxacin (MOX) and Gatifloxacin (GAT) on the Sensititre® Dried Susceptibility Plate Compared to NCCLS Microdilution Method

**Authors: T. Kelley<sup>1</sup>, S. Killian<sup>1</sup>, C.C. Knapp<sup>1</sup>, N. Holliday<sup>1</sup>, M.R. Jacobs<sup>2</sup>, A. Windau<sup>2</sup>, S. Bajaksouzian<sup>2</sup>, E. Smith Moland<sup>3</sup>, T. J. Lockhart<sup>3</sup>, L. B. Olson<sup>3</sup>, J. A. Black<sup>3</sup>**  
 1. TREK Diagnostic Systems, Cleveland, OH 2. University Hospital, Cleveland, OH 3. Creighton University, Omaha, NE



982 Keynote Circle, Suite 6  
 Cleveland, OH 44131  
 800.871.8909 • www.trekds.com

### ABSTRACT

**Background:** A multi-site evaluation was undertaken to compare the performance of MOX (Bayer, West Haven, CT) and GAT (BMSQ, Princeton, NJ) on the Sensititre 18 – 24 hour susceptibility dried plate to the NCCLS M7-broth microdilution reference method, testing both automated and manual reading methods. The range tested for both MOX and GAT was (0.004 – 16 µg/ml).

**Methods:** The clinical evaluation was conducted at three sites and consisted of 306 fresh clinical gram-negative isolates and 75 CDC challenge isolates, which included both *Enterobacteriaceae* and *non-Enterobacteriaceae*. Recommended quality control organisms were tested daily and were within the NCCLS expected quality control range.

**Results\*:** For MOX, when the automated reading method was compared to the reference method the total isolates' Essential Agreement (EA) and Categorical Agreement (CA) were 96.3% and 96.9%, respectively. When the manual reading method was compared to the reference method the EA and CA were 97.6% and 97.9%, respectively.

For GAT, when the automated reading method was compared to the reference method the total isolates' EA and CA were 97.4% and 98.4%, respectively. When the manual reading method was compared to the reference method the EA and CA were 97.6% and 99%, respectively.

For MOX, the reproducibility EA and CA were 97.3% and 92.0% for the automated method and 100% and 98.7% for the manual method. For GAT, the reproducibility EA and CA were 98.7% and 96% for the automated method and 100% and 97.3% for the manual method.

**Conclusions:** This multi-site evaluation indicates that the performance of MOX and GAT on the Sensititre 18-24 hour susceptibility system, using either the automated or manual reading method, was equivalent to its performance using the NCCLS M7-broth microdilution reference method.

\*Abstract has been modified

### PURPOSE OF THE STUDY

To evaluate the performance of Moxifloxacin and Gatifloxacin on the Sensititre 18 – 24 hour susceptibility plate compared to the NCCLS microdilution method (M7 – A6).

### MATERIALS & METHODS

**Organisms:** The testing at 3 sites consisted of the following:

- 306 clinical gram-negative isolates
- 75 CDC challenge isolates
- 25 reproducibility isolates
- Quality Control strains

#### Antibiotics

Antimicrobials Tested	Range Tested	Supplied By
Moxifloxacin (MOX)	0.004 – 16	Bayer
Gatifloxacin (GAT)	0.004 – 16	BMSQ

#### Susceptibility Testing Methods:

- Each isolate was tested using a Sensititre 18 – 24 hour susceptibility plates. The plates were set-up and tested according to the manufacturers' instructions.
- The reference plates were tested according to the microdilution methods published by the National Committee for Clinical Laboratory Standards (NCCLS, M7-A6).
- The approved primary "Indications for Use", and clinical significance of Gatifloxacin is for: *Escherichia coli*, *Klebsiella pneumoniae*, and *Proteus mirabilis*. *In vitro* data, without clinical correlation is provided for: *Acinetobacter lwoffii*, *Citrobacter koseri*, *Citrobacter freundii*, *Enterobacter cloacae*, *Klebsiella oxytoca*, *Morganella morganii*, and *Proteus vulgaris*.
- The approved primary "Indications for Use", and clinical significance of Moxifloxacin is for: *Klebsiella pneumoniae*. *In vitro* data, without clinical correlation is provided for: *Citrobacter freundii*, *Enterobacter cloacae*, *Klebsiella oxytoca*, and *Proteus mirabilis*.

Organisms	Isolates Tested		
	Clinical	Challenge	Total
<i>E. coli</i>	60	13	73
<i>Klebsiella species</i>	40	10	50
<i>Enterobacter species</i>	38	10	48
<i>Proteus species</i>	30	6	36
<i>Providencia species</i>	11	5	16
<i>Pseudomonas aeruginosa</i>	52	11	63
<i>Acinetobacter species</i>	15	6	21
*All others tested	60	14	74
<b>Total</b>	<b>306</b>	<b>75</b>	<b>381</b>

\**Citrobacter spp.*, *Morganella spp.*, *Serratia spp.*, *Aeromonas spp.*

### RESULTS

Gatifloxacin – Automated Read Method  
 MIC Distribution for Clinical and Challenge Isolates vs. Gatifloxacin

Test Results	Reference Results													
	0.004	0.008	0.015	0.03	0.06	0.12	0.25	0.5	1	2	4	8	16	>16
0.004														
0.008		4	2											
0.015			44	13	1									
0.03			5	55	20									
0.06			7	19	10									
0.12				4	5	35	10							
0.25				1	2	9	20	3						
0.5							2	16	10					
1							1	1	1	13	7			
2											13	6		
4												10		
8													12	5
16														4
>16														2
<b>TOTALS</b>	<b>0</b>	<b>4</b>	<b>51</b>	<b>88</b>	<b>47</b>	<b>55</b>	<b>33</b>	<b>20</b>	<b>23</b>	<b>20</b>	<b>16</b>	<b>12</b>	<b>11</b>	<b>9</b>

Categorical Interpretations: Susceptible ≤ 2, Intermediate = 4, Resistant ≥ 8

Gatifloxacin – Manual Read Method  
 MIC Distribution for Clinical and Challenge Isolates vs. Gatifloxacin

Test Results	Reference Results													
	0.004	0.008	0.015	0.03	0.06	0.12	0.25	0.5	1	2	4	8	16	>16
0.004														
0.008		4	2											
0.015			43	8	1									
0.03			6	62	15									
0.06			5	27	10									
0.12			4	4	39	9								
0.25			1		5	21	3	1						
0.5						2	15	10						
1							1	2	11	9				
2								1	11	3				
4											12			
8												1	10	6
16													2	4
>16														1
<b>TOTALS</b>	<b>0</b>	<b>4</b>	<b>51</b>	<b>80</b>	<b>47</b>	<b>55</b>	<b>33</b>	<b>20</b>	<b>23</b>	<b>20</b>	<b>16</b>	<b>12</b>	<b>11</b>	<b>9</b>

Categorical Interpretations: Susceptible ≤ 2, Intermediate = 4, Resistant ≥ 8

Reproducibility Isolates vs. Gatifloxacin / Moxifloxacin  
 Reference Agreement within ±1 Dilution of Modal MIC

Site	# Agreement		Total Isolates Tested	Manual %	Autoread %
	Manual	Autoread			
TREK	25/25	25/25	25/25	100/100	100/100
Creighton	25/25	25/25	25/25	100/100	100/100
UH	25/25	24/23	25/25	100/100	96.0/92.0
Combined	75/75	74/73	75/75	100/100	98.7/97.3

Moxifloxacin – Automated Read Method  
 MIC Distribution for Clinical and Challenge Isolates vs. Moxifloxacin

Test Results	Reference Results													
	0.004	0.008	0.015	0.03	0.06	0.12	0.25	0.5	1	2	4	8	16	>16
0.004														
0.008														
0.015				5	5	1								
0.03				3	42	14								
0.06				1	9	59	14							
0.12					5	16	9	1						
0.25						6	8	35	2	1				
0.5							7	26	9					
1								1		2	14	13		
2									1			13	6	
4										1			7	4
8													1	9
16														11
>16														9
<b>TOTALS</b>	<b>0</b>	<b>0</b>	<b>9</b>	<b>56</b>	<b>85</b>	<b>40</b>	<b>52</b>	<b>31</b>	<b>24</b>	<b>26</b>	<b>14</b>	<b>13</b>	<b>20</b>	<b>11</b>

Categorical Interpretations: Susceptible ≤ 2, Intermediate = 4, Resistant ≥ 8

Moxifloxacin – Manual Read Method  
 MIC Distribution for Clinical and Challenge Isolates vs. Moxifloxacin

Test Results	Reference Results													
	0.004	0.008	0.015	0.03	0.06	0.12	0.25	0.5	1	2	4	8	16	>16
0.004														
0.008														
0.015				5	4									
0.03				3	45	9								
0.06				1	7	65	10							
0.12						5	22	7						
0.25						6	7	35	2					
0.5							9	26	7					
1								1		3	16	8		
2											1	18	3	
4												10	3	
8													1	10
16														4
>16														1
<b>TOTALS</b>	<b>0</b>	<b>0</b>	<b>9</b>	<b>56</b>	<b>85</b>	<b>40</b>	<b>52</b>	<b>31</b>	<b>24</b>	<b>26</b>	<b>14</b>	<b>13</b>	<b>15</b>	<b>16</b>

Categorical Interpretations: Susceptible ≤ 2, Intermediate = 4, Resistant ≥ 8

Gatifloxacin versus Clinical and Challenge Isolates

Read Method	Total Isolates Tested	% Essential Agreement	% Categorical Agreement	Number Resistant	Very Major Errors	Major Errors	Minor Errors
Manual	381	97.6	99.0	32	0	0	4
Automated	381	97.4	98.4	32	0	0	6

Moxifloxacin versus Clinical and Challenge Isolates

Read Method	Total Isolates Tested	% Essential Agreement	% Categorical Agreement	Number Resistant	Very Major Errors	Major Errors	Minor Errors
Manual	381	97.6	97.9	44	0	0	8
Automated	381	96.3	96.9	44	0	0	12

### CONCLUSIONS

This investigation compared the 18 – 24 hour Sensititre® susceptibility plate with the NCCLS reference microdilution frozen plate (M7 – A6). The Sensititre plate demonstrated a high level of agreement and was very reproducible.

#### Clinical Isolates

- Gatifloxacin:**
- The overall essential agreement for Gatifloxacin, within a +/- one-well dilution range, was 97.1% for the manual method and 96.7% for the autoread method.
  - The overall categorical agreement for Gatifloxacin was 99.0% for the manual method and 99.0% for the autoread method.

- Moxifloxacin:**
- The overall essential agreement for Moxifloxacin, within a +/- one-well dilution range, was 97.1% for the manual method and 95.8% for the autoread method.
  - The overall categorical agreement for Moxifloxacin was 98.4% for the manual method and 97.4% for the autoread method.

#### CDC Challenge Organisms

- Gatifloxacin:**
- The overall essential agreement for Gatifloxacin, within a +/- one-well dilution range, was 100% for the manual method and 100% for the autoread method.
  - The overall categorical agreement for Gatifloxacin was 98.7% for the manual method and 96.0% for the autoread method.

- Moxifloxacin:**
- The overall essential agreement for Moxifloxacin, within a +/- one-well dilution range, was 100% for the manual method and 98.7% for the autoread method. The overall categorical agreement for Moxifloxacin was 96.0% for the manual method and 94.7% for the autoread method.

#### Reproducibility

- Gatifloxacin:**
- For Gatifloxacin the Interlaboratory essential agreement reproducibility was 100% for the manual read method and 98.7% for the autoread method.
  - The Interlaboratory categorical agreement reproducibility was 97.3% for the manual read method and 96.0% for the autoread method.

- Moxifloxacin:**
- For Moxifloxacin the Interlaboratory essential agreement reproducibility was 100% for the manual read method and 97.3% for the autoread method.
  - The Interlaboratory categorical agreement reproducibility was 98.7% for the manual read method and 92.0% for the autoread method.