

# Reproducibility of mycobacterial QC strains in the Sensititre® rapid growing and slowly growing nontuberculosis mycobacteria broth microdilution susceptibility procedure

Brimecombe, M., P. Stagg, G. Hammerton, S. Killian, C. Knapp and R. Grist  
TREK Diagnostics Systems, East Grinstead, UK and Cleveland Ohio, USA

Trek Diagnostic Inc.  
982 Keynote Circle, Suite 6  
Cleveland, OH 44131  
800.871.8909 • www.trekds.com

## ABSTRACT

### OBJECTIVE

Presently there are no NCCLS QC ranges for *M. avium* and only a few for rapid growing mycobacteria. The purpose of this study was to measure reproducibility with *M. peregrinum* ATCC 700686, *M. smegmatis* ATCC 19420 and *M. avium* ATCC 700898 on commercial Sensititre dried susceptibility panels. For rapid growers, the effect of incubation time and the addition of saline / 0.2% tween with or without glass beads were also investigated. The antimicrobials evaluated were: Amikacin, Ofloxacin, Ansamycin, Streptomycin, Clofazimine, Kanamycin, Capreomycin, Ofloxacin, Rifampicin, Ethambutol and Isoniazid.

### METHOD; rapid growers

A 0.5 McFarland standard was prepared from colonies on a 48hr old TSA / blood plate. Suspensions were vortexed for 5 minutes allowing large clumps to settle. 50 µL were transferred to 10 mL Sensititre Mueller Hinton broth. 100 µL were dosed into each well in the plate. Plates were sealed and incubated at 30°C. Plates were examined for visual growth after 72, 96 and 120 hours. Two sites tested 25 and 50 replicate MICs for each isolate. An additional 10 results were collected for each of the three suspending media.

### METHOD; slow growers

The set up procedure was the same except that 7H9 broth or Sensititre Mueller Hinton broth were supplemented with OADC. Plates were sealed and incubated at 35°C and read at 10 to 14 days depending upon extent of growth. A total of 50 replicates MICs were collected.

### RESULTS; rapid growers

All onscale MICs fell within one doubling dilution of the mode. The maximum shift in MIC over time was one dilution. The presence of tween resulted in a one well shift in MIC of Ansmycin and Clofazimine with ATCC 19420 and Clofazimine with ATCC 700686.

### RESULTS; slow growers

All MICs fell within 4 dilutions with the exception of Ofloxacin. Weak growth requiring 14 days incubation tended overall to give lower MICs but fell within the 4 well range.

### CONCLUSION

MICs were highly reproducible allowing tentative QC ranges to be set for testing Sensititre plates. MICs were not affected by the length of incubation. Saline / tween showed notable effects on the MICs with two of the drugs with rapid growers.

## INTRODUCTION

*Mycobacterium avium* Complex (MAC) causes a wide range of infections including disseminated disease in patients with AIDS and invasive pulmonary disease. The rapidly growing pathogenic mycobacteria cause several forms of disease of varying severity, most commonly skin and soft tissue infections but also skeletal, pulmonary and disseminated disease. Antimicrobial susceptibility testing of clinically significant isolates is recommended. NCCLS (1) has published a recommended procedure but there is currently only a limited number of drug QC ranges for testing rapidly growing pathogenic mycobacteria and MAC. Other antimicrobials have to be quality controlled with non-fastidious organisms listed in Table 3 of NCCLS M100-Performance Standards for Antimicrobial Susceptibility testing (2).

The interlaboratory reproducibility of Sensititre dried microdilution panels for susceptibility testing of rapidly growing mycobacteria (3) and *Mycobacterium avium* Complex (4) has been evaluated and shown to be acceptable. The purpose of this study was to measure reproducibility with rapid grower QC strains, *M. peregrinum* ATCC 700686, *M. smegmatis* ATCC 19420 and MAC QC strain *M. avium* ATCC 700898 on Sensititre dried susceptibility panels. For rapid growers, the effect of incubation time and the addition of saline / 0.2% Tween with or without glass beads to suspending media was also investigated. MAC testing compared MICs using Mueller Hinton broth with OADC compared to 7HSF broth.

## METHODS

### Rapid growing Mycobacteria

0.5 McFarland standard prepared from colonies on a 48hr old TSA / blood plate.

Suspensions vortexed in sterile water for 5 minutes allowing large clumps to settle

50 µL transferred to 10 mL Sensititre Mueller Hinton broth To obtain a final organism concentration of 1 X 10<sup>8</sup> CFU/mL to 5 X 10<sup>8</sup> CFU/mL

100 µL dosed into each well of the plate. Plates sealed and incubated at 30°C

Plates examined for visual growth after 72 hours.15 plates had incubation extended to 96 and 120 hours. Two sites collected 25 and 50 replicate results. Testing was performed over 5 days. An additional 10 results were collected from suspensions vortexed in saline containing 0.2% Tween 80 with or without glass beads.

Plates inoculated with suspensions prepared in Saline/Tween with or without beads were incubated for 72 hours

### M. avium

The set up procedure same as for rapid growers except for use of 7H9 broth supplemented with 5% OADC (7HSF) or Sensititre Mueller Hinton broth pH 7.3 supplemented with 5% OADC (MHB-OADC)

Plates inoculated with MHB-OADC broth were sealed and incubated at 35°C and read at 10 or 14 days depending upon extent of growth. A total of 50 replicates from one site MICs were collected.

Ten plates dosed with 7HSF broth were read at 6 to 21 days.

## RESULTS

Table 1. Effect of incubation time and suspension medium on reproducibility of *M. peregrinum* ATCC 700686 MICs

Antibiotic/Experimental variable	Number of results at MIC (µg/mL)
Amikacin	0.25 0.5 1 2 4
72h	66 9
96h	5 10
120h	1 14
Saline/Tween	10
Saline/Tween beads	9 1
Ansamycin	0.5 1 2 4 8
72h	1 46 28
96h	4 4 11
120h	15 15
Saline/Tween	6 4
Saline/Tween beads	6 4
Capreomycin	1.25 2.5 5 10 20
72h	13 53 9
96h	7 8 7
120h	5 15
Saline/Tween	5 5
Saline/Tween beads	4 4 6
Clofazimine	0.12 0.25 0.5 1 2
72h	45 25 5
96h	7 3
120h	1 12 2
Saline/Tween	10
Saline/Tween beads	8 2
Ethambutol	2 4 8 16 32
72h	2 16 54 3
96h	10 5
120h	4 11
Saline/Tween	10
Saline/Tween beads	10
Ethionamide	>20
72h	75
96h	15 15
120h	10
Saline/Tween	10
Saline/Tween beads 10	10
Isoniazid	<4
72h	75
96h	15 34 24 16
120h	10 5
Saline/Tween	10
Saline/Tween beads 10	10
Kanamycin	1.25 2.5 5 10 20
72h	7 46 22 3
96h	14 1
120h	14 1
Saline/Tween	10
Saline/Tween beads	10
Ofloxacin	<0.25 0.5 1 2 4
72h	85 10
96h	14 1
120h	14 1
Saline/Tween	10
Saline/Tween beads	10
Rifampicin	>16
72h	75
96h	15 10
120h	10
Saline/Tween	10
Saline/Tween beads 10	10
Streptomycin	16 32 >32
72h	1 40 34
96h	1 14
120h	1 15
Saline/Tween	4 6
Saline/Tween beads	2 8

Table 2. Effect of incubation time and suspension medium on reproducibility of *M. smegmatis* ATCC 19420 MICs

Antibiotic/Experimental variable	Number of results at MIC (µg/mL)
Amikacin	0.12 0.25 0.5 1 2
72h	72 3
96h	15
120h	15
Saline/Tween	10
Saline/Tween beads	10
Ansamycin	0.5 1 2 4 8
72h	2 72 1
96h	14 11
120h	15 13 2
Saline/Tween	3 2 5
Saline/Tween beads	1 5 4
Capreomycin	0.06 1.25 2.5 5 10
72h	15 60
96h	1 14
120h	1 15
Saline/Tween	4 6
Saline/Tween beads	3 7
Clofazimine	0.06 0.12 0.25 0.5 1
72h	22 47 6
96h	2 12 1
120h	1 12 2
Saline/Tween	2 8
Saline/Tween beads	10
Ethambutol	<0.5
72h	75
96h	15 14
120h	15
Saline/Tween	10
Saline/Tween beads	10
Ethionamide	20 >20 40 80 160
72h	16 69 11* 44*
96h	15
120h	10
Saline/Tween	2 7
Saline/Tween beads	5 5
*Note: Different panel with extended concentration range to bring MICs onscale	
Isoniazid	0.5 1 2 4 >4
72h	75
96h	1 34 24 16
120h	10 5
Saline/Tween	2 3
Saline/Tween beads	4 1
Kanamycin	<0.06 1.25 2.5 1 2
72h	75
96h	15 14
120h	15
Saline/Tween	10
Saline/Tween beads	10
Ofloxacin	0.25 0.5 1 2 4
72h	71 4
96h	15
120h	15
Saline/Tween	10
Saline/Tween beads	10
Rifampicin	>16
72h	75
96h	15 10
120h	10
Saline/Tween	10
Saline/Tween beads	10
Streptomycin	0.25 0.5 1 2 4
72h	28 46
96h	6 9
120h	4 11
Saline/Tween	9 1
Saline/Tween beads	10 10

Table 3. Effect of incubation time and broth type on reproducibility of *M. avium* ATCC 700898 MICs

Broth	Incubation time (days)	Number of Amikacin MICs (µg/mL)				
MHB OADC	10	1	2	4	8	>8
MHB OADC	14	3	21	5	1	
7HSF	6			4	6	
7HSF	10			4	5	1
7HSF	14			4	4	1
7HSF	21			2	3	1
Broth	Incubation time (days)	Number of Ansamycin MICs (µg/mL)				
MHB OADC	10	0.06	0.12	0.25	0.5	1 2
MHB OADC	14			12	8	
7HSF	6			10	18	6 2
7HSF	7			4	20	
7HSF	10			1	4	1
7HSF	14			1	3	1
7HSF	21			2	3	
Broth	Incubation time (days)	Number of Capreomycin MICs (µg/mL)				
MHB OADC	10	5	10	20		>20
MHB OADC	14	5	15	9		
7HSF	6			1	9	
7HSF	10			1	4	1
7HSF	14			3	5	
7HSF	21			2	5	
Broth	Incubation time (days)	Number of Clofazimine MICs (µg/mL)				
MHB OADC	10	≤ 0.06	0.12	0.25	0.5	1 2
MHB OADC	14		11			
7HSF	6			1	9	
7HSF	7				5	
7HSF	10				5	2 3
7HSF	14				2	4
7HSF	21				3	2
Broth	Incubation time (days)	Number of Ethambutol MICs (µg/mL)				
MHB OADC	10	4	8	16	32	64 >64
MHB OADC	14		9			
7HSF	6		26			
7HSF	7					10
7HSF	10			4	1	2 6
7HSF	14			4	4	1 3
7HSF	21			3	3	2
Broth	Incubation time (days)	Number of Ethionamide MICs (µg/mL)				
MHB OADC	10	1.25	2.5	5	10	20 >20
MHB OADC	14		11			
7HSF	6		12	6	1	
7HSF	7		2	12	15	1
7HSF	10		5	1	1	2 2 6
7HSF	14			2	2	2 1
7HSF	21			2	2	5 3
Broth	Incubation time (days)	Number of Isoniazid MICs (µg/mL)				
MHB OADC	10	1	2	4	>4	
MHB OADC	14				20	
7HSF	6				30	
7HSF	7				10	
7HSF	10				19	
7HSF	14				2	2
7HSF	21				5	4 1
Broth	Incubation time (days)	Number of Kanamycin MICs (µg/mL)				
MHB OADC	10	1.25	2.5	5	10	20
MHB OADC	14			1	19	2
7HSF	6		3	6	7	3
7HSF	7				3	2
7HSF	10				3	2
7HSF	14				4	1
7HSF	21				5	
Broth	Incubation time (days)	Number of Ofloxacin MICs (µg/mL)				
MHB OADC	10	2	4	8	16	>16
MHB OADC	14		8	18	4	
7HSF	6				9	
7HSF	7			2	3	
7HSF	10				4	1
7HSF	14				4	1
7HSF	21				3	2
Broth	Incubation time (days)	Number of Rifampicin MICs (µg/mL)				
MHB OADC	10	1	2	4	8	16 >16
MHB OADC	14				12	
7HSF	6				4	8
7HSF	7			2	8	
7HSF	10			4	1	
7HSF	14			1	3	1
7HSF	21			2	2	1
Broth	Incubation time (days)	Number of Streptomycin MICs (µg/mL)				
MHB OADC	10	2	4	8	16	32
MHB OADC	14		4	8	12	6
7HSF	6			7	3	
7HSF	7			1	4	1
7HSF	10			1	4	1
7HSF	14			1	4	1
7HSF	21			1	4	1

## CONCLUSIONS

### Rapid growers

Tables 1 and 2 show excellent reproducibility between sites

Onscale MICs fell within one doubling dilution of the mode with exception of Ethambutol with *M. peregrinum* ATCC 700686 (1.9% outliers) and Isoniazid with *M. smegmatis* ATCC 19420 (0.9% outliers)

Use of Saline and Tween to prepare inoculum resulted in a one well shift in the MIC for Anamycin and Clofazimine with *M. smegmatis* ATCC 19420 and Clofazimine with *M. peregrinum* ATCC 700686 on panels read after 72 hours incubation

Extending incubation to from 72 hours to 120 hours increased the modal MIC for all drugs giving onscale MICs by one dilution with *M. peregrinum* ATCC 700686. No shift was observed with *M. smegmatis* ATCC 19420 apart from Isoniazid.

Table 4 shows tentative QC ranges based upon combining data from Tables 1 and 2 with other published information.

### Slow Growers

Onscale MICs fell within 4 dilutions with the exception of Ofloxacin with MHB-OADC broth when incubated for 10 to 14 days

Weak growth requiring 14 days incubation gave lower MICs but fell with the 4 well range for Amikacin, Capreomycin, Clofazimine, Ethambutol, Ofloxacin and Streptomycin

*M. avium* ATCC 700898 with 7HSF broth gave poorly reproducible results with Ethambutol, Clofazimine and Ethionamide