

**SENSITITRE™ GRAM NEGATIVE PLATE FORMAT**

Plate Code: **GN4F**

	1	2	3	4	5	6	7	8	9	10	11	12
<b>A</b>	AMI 8	TGC 1	LEVO 1	DOR 0.5	ETP 0.25	IMI 0.5	MERO 0.5	FAZ 1	TAZ 1	AZT 1	FEP 4	AXO 0.5
<b>B</b>	AMI 16	TGC 2	LEVO 2	DOR 1	ETP 0.5	IMI 1	MERO 1	FAZ 2	TAZ 2	AZT 2	FEP 8	AXO 1
<b>C</b>	AMI 32	TGC 4	LEVO 4	DOR 2	ETP 1	IMI 2	MERO 2	FAZ 4	TAZ 4	AZT 4	FEP 16	AXO 2
<b>D</b>	P/T4 8/4	TGC 8	LEVO 8	DOR 4	ETP 2	IMI 4	MERO 4	FAZ 8	TAZ 8	AZT 8	FEP 32	AXO 4
<b>E</b>	P/T4 16/4	TIM2 8/2	NIT 32	MIN 1	ETP 4	IMI 8	MERO 8	FAZ 16	TAZ 16	AZT 16	CIP 0.5	AXO 8
<b>F</b>	P/T4 32/4	TIM2 16/2	NIT 64	MIN 2	ETP 8	PIP 16	GEN 2	TOB 2	A/S2 4/2	AMP 8	CIP 1	AXO 16
<b>G</b>	P/T4 64/4	TIM2 32/2	TET 4	MIN 4	SXT 2/38	PIP 32	GEN 4	TOB 4	A/S2 8/4	AMP 16	CIP 2	AXO 32
<b>H</b>	P/T4 128/4	TIM2 64/2	TET 8	MIN 8	SXT 4/76	PIP 64	GEN 8	TOB 8	A/S2 16/8	POS	POS	POS

**ANTIMICROBICS**

AMI	Amikacin
P/T4	Piperacillin / tazobactam constant 4
TGC	Tigecycline
TIM2	Ticarcillin / clavulanic acid constant 2
LEVO	Levofloxacin
NIT	Nitrofurantoin
TET	Tetracycline
DOR	Doripenem
MIN	Minocycline
ETP	Ertapenem
SXT	Trimethoprim / sulfamethoxazole
IMI	Imipenem
PIP	Piperacillin
MERO	Meropenem
GEN	Gentamicin
FAZ	Cefazolin
TOB	Tobramycin
TAZ	Ceftazidime
A/S2	Ampicillin / sulbactam 2:1 ratio
AZT	Aztreonam
AMP	Ampicillin
POS	Positive Control
FEP	Cefepime
CIP	Ciprofloxacin
AXO	Ceftriaxone