



## SENSITITRE™ GRAM NEGATIVE URINE PLATE FORMAT

Plate Code: **GNUR2F**

	1	2	3	4	5	6	7	8	9	10	11	12
<b>A</b>	NEG	POS	POS	POS	NEG	POS	POS	POS	NEG	POS	POS	POS
<b>B</b>	SXT 0.5/9.5	SXT 2/38	SXT 4/76	FEP 32	SXT 0.5/9.5	SXT 2/38	SXT 4/76	FEP 32	SXT 0.5/9.5	SXT 2/38	SXT 4/76	FEP 32
<b>C</b>	NIT 32	NIT 64	NIT 128	FEP 16	NIT 32	NIT 64	NIT 128	FEP 16	NIT 32	NIT 64	NIT 128	FEP 16
<b>D</b>	CIP 1	CIP 2	CIP 4	FEP 8	CIP 1	CIP 2	CIP 4	FEP 8	CIP 1	CIP 2	CIP 4	FEP 8
<b>E</b>	AUG2 8/4	AUG2 16/8	AUG2 32/16	FEP 4	AUG2 8/4	AUG2 16/8	AUG2 32/16	FEP 4	AUG2 8/4	AUG2 16/8	AUG2 32/16	FEP 4
<b>F</b>	AMP 8	AMP 16	AMP 32	GEN 8	AMP 8	AMP 16	AMP 32	GEN 8	AMP 8	AMP 16	AMP 32	GEN 8
<b>G</b>	CAR 16	CAR 32	CAR 64	GEN 4	CAR 16	CAR 32	CAR 64	GEN 4	CAR 16	CAR 32	CAR 64	GEN 4
<b>H</b>	AXO 8	AXO 16	AXO 32	AXO 64	AXO 8	AXO 16	AXO 32	AXO 64	AXO 8	AXO 16	AXO 32	AXO 64

### ANTIMICROBICS

NEG	Negative Control
SXT	Trimethoprim / sulfamethoxazole
NIT	Nitrofurantoin
CIP	Ciprofloxacin
AUG2	Amoxicillin / clavulanic acid 2:1 ratio
AMP	Ampicillin
CAR	Carbenicillin
AXO	Ceftriaxone
POS	Positive Control
FEP	Cefepime
GEN	Gentamicin