

SENSITITRE™ GRAM POSITIVE PLATE FORMAT

Plate Code: **GPALL2F**

	1	2	3	4	5	6	7	8	9	10	11	12
A	CPT 0.12	CPT 0.25	CPT 0.5	CPT 1	CPT 2	CPT 4	AMP 0.25	AMP 0.5	AMP 1	AMP 2	AMP 4	AMP 8
B	CHL 2	CHL 4	CHL 8	CHL 16	AXO 8	AXO 16	AXO 32	AXO 64	DAP 0.5	DAP 1	DAP 2	DAP 4
C	CLI 0.5	CLI 1	CLI 2	GEN 4	GEN 8	LEVO 0.25	LEVO 0.5	LEVO 1	LEVO 2	LEVO 4	DT1	DT2
D	LZD 1	LZD 2	LZD 4	LZD 8	VAN 0.25	VAN 0.5	VAN 1	VAN 2	VAN 4	VAN 8	VAN 16	FOXS 6
E	PEN 0.12	PEN 0.25	PEN 0.5	PEN 1	PEN 2	PEN 4	PEN 8	TLV 0.06	TLV 0.12	TLV 0.25	TLV 0.5	TLV 1
F	MXF 0.25	MXF 0.5	MXF 1	MXF 2	MXF 4	OXA+ 0.25	OXA+ 0.5	OXA+ 1	OXA+ 2	SXT 2/38	SXT 4/76	POS
G	RIF 1	RIF 2	TGC 0.03	TGC 0.06	TGC 0.12	TGC 0.25	TGC 0.5	NIT 32	NIT 64	GEN 500	STR 1000	POS
H	TET 2	TET 4	TET 8	ERY 0.5	ERY 1	ERY 2	ERY 4	CIP 1	CIP 2	CIP 4	NEG	POS

ANTIMICROBICS

CHL	Chloramphenicol
DAP	Daptomycin
AMP	Ampicillin
PEN	Penicillin
LZD	Linezolid
MXF	Moxifloxacin
RIF	Rifampin
TET	Tetracycline
TGC	Tigecycline
ERY	Erythromycin
CPT	Ceftaroline
LEVO	Levofloxacin
VAN	Vancomycin
OXA+	Oxacillin+2%NaCl
AXO	Ceftriaxone
TLV	Telavancin
NIT	Nitrofurantoin
CIP	Ciprofloxacin
CLI	Cilindamycin
GEN	Gentamicin
SXT	Trimethoprim / sulfamethoxazole
DT1	D Test 1
STR	Streptomycin
NEG	Negative Control
FOXS	Cefoxitin screen
DT2	D Test 2
POS	Positive Control