

SENSITITRE™ GRAM POSITIVE PLATE FORMAT

Plate Code: **GPN3F**

	1	2	3	4	5	6	7	8	9	10	11	12
A	ERY 0.25	ERY 0.5	ERY 1	ERY 2	ERY 4	CLI 0.12	CLI 0.25	CLI 0.5	CLI 1	CLI 2	GEN 500	STR 1000
B	SYN 0.12	SYN 0.25	SYN 0.5	SYN 1	SYN 2	SYN 4	DAP 0.25	DAP 0.5	DAP 1	DAP 2	DAP 4	DAP 8
C	VAN 1	VAN 2	VAN 4	VAN 8	VAN 16	VAN 32	VAN 64	VAN 128	TET 2	TET 4	TET 8	TET 16
D	AMP 0.12	AMP 0.25	AMP 0.5	AMP 1	AMP 2	AMP 4	AMP 8	AMP 16	GEN 2	GEN 4	GEN 8	GEN 16
E	RIF 0.5	LEVO 0.25	LEVO 0.5	LEVO 1	LEVO 2	LEVO 4	LEVO 8	LZD 0.5	LZD 1	LZD 2	LZD 4	LZD 8
F	RIF 1	PEN 0.06	PEN 0.12	PEN 0.25	PEN 0.5	PEN 1	PEN 2	PEN 4	PEN 8	CIP 0.5	CIP 1	CIP 2
G	RIF 2	SXT 1/19	SXT 2/38	SXT 4/76	AXO 8	AXO 16	AXO 32	AXO 64	GAT 1	GAT 2	GAT 4	GAT 8
H	RIF 4	SXT 0.5/9.5	OXA+ 0.25	OXA+ 0.5	OXA+ 1	OXA+ 2	OXA+ 4	OXA+ 8	NEG	POS	POS	POS

ANTIMICROBICS

ERY	Erythromycin
CLI	Clindamycin
SYN	Quinupristin / dalfopristin
DAP	Daptomycin
VAN	Vancomycin
TET	Tetracycline
AMP	Ampicillin
GEN	Gentamicin
LEVO	Levofloxacin
LZD	Linezolid
AXO	Ceftriaxone
STR	Streptomycin
PEN	Penicillin
RIF	Rifampin
GAT	Gatifloxacin
CIP	Ciprofloxacin
SXT	Trimethoprim / sulfamethoxazole
OXA+	Oxacillin+2%NaCL
POS	Positive Control
NEG	Negative Control