

SENSITITRE™ STREP PLATE FORMAT

Plate Code: **STP6F**

	1	2	3	4	5	6	7	8	9	10	11	12
A	MXF 1	MXF 2	MXF 4	MXF 8	PEN 0.03	PEN 0.06	PEN 0.12	PEN 0.25	PEN 0.5	PEN 1	PEN 2	PEN 4
B	LEVO 0.5	LEVO 1	LEVO 2	LEVO 4	MERO 0.25	MERO 0.5	MERO 1	MERO 2	AZI 0.25	AZI 0.5	AZI 1	AZI 2
C	TET 1	TET 2	TET 4	TET 8	ETP 0.5	ETP 1	ETP 2	ETP 4	ERY 0.25	ERY 0.5	ERY 1	ERY 2
D	FUR 0.5	FUR 1	FUR 2	FUR 4	AUG2 2/1	AUG2 4/2	AUG2 8/4	AUG2 16/8	SXT 0.5/9.5	SXT 1/19	SXT 2/38	SXT 4/76
E	AXO 0.12	AXO 0.25	AXO 0.5	AXO 1	AXO 2	LZD 0.25	LZD 0.5	LZD 1	LZD 2	LZD 4	VAN 0.5	VAN 1
F	FOT 0.12	FOT 0.25	FOT 0.5	FOT 1	FOT 2	FOT 4	CLI 0.12	CLI 0.25	CLI 0.5	CLI 1	VAN 2	VAN 4
G	DAP 0.06	DAP 0.12	DAP 0.25	DAP 0.5	DAP 1	DAP 2	FEP 0.5	FEP 1	FEP 2	FEP 4	FEP 8	POS
H	CHL 1	CHL 2	CHL 4	CHL 8	CHL 16	CHL 32	TGC 0.015	TGC 0.03	TGC 0.06	TGC 0.12	POS	POS

ANTIMICROBICS

MXF	Moxifloxacin
LEVO	Levofloxacin
TET	Tetracycline
FUR	Cefuroxime
AXO	Ceftriaxone
FOT	Cefotaxime
DAP	Daptomycin
CHL	Chloramphenicol
PEN	Penicillin
MERO	Meropenem
ETP	Ertapenem
AUG2	Amoxicillin / clavulanic acid 2:1 ratio
LZD	Linezolid
CLI	Clindamycin
FEP	Cefepime
TGC	Tigecycline
AZI	Azithromycin
ERY	Erythromycin
SXT	Trimethoprim / sulfamethoxazole
VAN	Vancomycin
POS	Positive Control