



SENSITITRE AQUATIC PLATE FORMAT

Plate Code: **AQUATIC**

Plate Type: **MIC**

	1	2	3	4	5	6	7	8	9	10	11	12
A	ENRO 1	ENRO 0.5	ENRO 0.25	ENRO 0.12	ENRO 0.06	ENRO 0.03	ENRO 0.015	ENRO 0.008	ENRO 0.004	ENRO 0.002	GEN 4	SXT 1/19
B	AMP 16	AMP 8	AMP 4	AMP 2	AMP 1	AMP 0.5	AMP 0.25	AMP 0.12	AMP 0.06	AMP 0.03	GEN 2	SXT 0.5/9.5
C	OXY 8	OXY 4	OXY 2	OXY 1	OXY 0.5	OXY 0.25	OXY 0.12	OXY 0.06	OXY 0.03	OXY 0.015	GEN 1	SXT 0.25/4.75
D	ERY 128	ERY 64	ERY 32	ERY 16	ERY 8	ERY 4	ERY 2	ERY 1	ERY 0.5	ERY 0.25	GEN 0.5	SXT 0.12/2.38
E	FFN 16	FFN 8	FFN 4	FFN 2	FFN 1	FFN 0.5	FFN 0.25	FFN 0.12	FFN 0.06	FFN 0.03	GEN 0.25	SXT 0.06/1.19
F	FLUQ 4	FLUQ 2	FLUQ 1	FLUQ 0.5	FLUQ 0.25	FLUQ 0.12	FLUQ 0.06	FLUQ 0.03	FLUQ 0.015	FLUQ 0.008	GEN 0.12	SXT 0.03/0.59
G	PRI 76/4	PRI 38/2	PRI 19/1	PRI 9.5/0.5	PRI 4.75/0.25	PRI 2.38/0.12	PRI 1.19/0.06	PRI 0.59/0.03	PRI 0.30/0.015	PRI 0.008/0.15	GEN 0.06	SXT 0.015/0.3
H	OXO 2	OXO 1	OXO 0.5	OXO 0.25	OXO 0.12	OXO 0.06	OXO 0.03	OXO 0.015	OXO 0.008	OXO 0.004	POS	NEG

ANTIMICROBICS

- SXT Trimethoprim / sulfamethoxazole
- GEN Gentamicin
- ENRO Enrofloxacin
- AMP Ampicillin
- OXY Oxytetracycline
- ERY Erythromycin
- FFN Florfenicol
- FLUQ Flumequine
- PRI Ormetoprim/sulphadimethoxine
- OXO Oxolinic Acid
- POS Positive Control
- NEG Negative Control ****NORMAL CONTROL**