



### SENSITITRE AVIAN PLATE FORMAT

Plate Code: **AVIAN1F**

Plate Type: **MIC**

	1	2	3	4	5	6	7	8	9	10	11	12
<b>A</b>	ENRO 2	ENRO 1	ENRO 0.5	ENRO 0.25	ENRO 0.12	SPE 64	SDM 256	FFN 8	PEN 8	STR 1024	NOV 4	CLI 4
<b>B</b>	GEN 8	GEN 4	GEN 2	GEN 1	GEN 0.5	SPE 32	SDM 128	FFN 4	PEN 4	STR 512	NOV 2	CLI 2
<b>C</b>	TIO 4	TIO 2	TIO 1	TIO 0.50	TIO 0.25	SPE 16	SDM 64	FFN 2	PEN 2	STR 256	NOV 1	CLI 1
<b>D</b>	NEO 32	NEO 16	NEO 8	NEO 4	NEO 2	SPE 8	SDM 32	FFN 1	PEN 1	STR 128	NOV 0.5	CLI 0.5
<b>E</b>	ERY 4	ERY 2	ERY 1	ERY 0.5	ERY 0.25	ERY 0.12	SXT 2/38	STZ 256	PEN 0.5	STR 64	TYLT 20	NEG
<b>F</b>	OXY 8	OXY 4	OXY 2	OXY 1	OXY 0.5	OXY 0.25	SXT 1/19	STZ 128	PEN 0.25	STR 32	TYLT 10	POS
<b>G</b>	TET 8	TET 4	TET 2	TET 1	TET 0.5	TET 0.25	SXT 0.5/9.5	STZ 64	PEN 0.12	STR 16	TYLT 5	POS
<b>H</b>	AMOX 16	AMOX 8	AMOX 4	AMOX 2	AMOX 1	AMOX 0.5	AMOX 0.25	STZ 32	PEN 0.06	STR 8	TYLT 2.5	POS

**ANTIMICROBICS**

ENRO	Enrofloxacin
GEN	Gentamicin
TIO	Ceftiofur
NEO	Neomycin
ERY	Erythromycin
OXY	Oxytetracycline
TET	Tetracycline
AMOX	Amoxicillin
SPE	Spectinomycin
SDM	Sulphadimethoxine
SXT	Trimethoprim / sulfamethoxazole
FFN	Florfenicol
STZ	Sulphathiazole
PEN	Penicillin
STR	Streptomycin
NOV	Novobiocin
TYLT	Tylosin tartrate
CLI	Clindamycin
POS	Postive Control
NEG	Negative Control