



SENSITITRE CUSTOM PLATE FORMAT

Plate Code: **BOPO6F**

	1	2	3	4	5	6	7	8	9	10	11	12
A	TIO 8	TIA 32	CTET 8	OXY 8	PEN 8	AMP 16	DANO 1	SXT 2/38	TYLT 4	TUL 4	CLI 16	SDM 256
B	TIO 4	TIA 16	CTET 4	OXY 4	PEN 4	AMP 8	DANO 0.5	SPE 64	TYLT 2	TUL 2	CLI 8	ENRO 2
C	TIO 2	TIA 8	CTET 2	OXY 2	PEN 2	AMP 4	DANO 0.25	SPE 32	TYLT 1	TUL 1	CLI 4	ENRO 1
D	TIO 1	TIA 4	CTET 1	OXY 1	PEN 1	AMP 2	DANO 0.12	SPE 16	TYLT 0.5	TIL 64	CLI 2	ENRO 0.5
E	TIO 0.5	TIA 2	CTET 0.5	OXY 0.5	PEN 0.5	AMP 1	NEO 32	SPE 8	TUL 64	TIL 32	CLI 1	ENRO 0.25
F	TIO 0.25	TIA 1	TIA 0.5	GEN 16	PEN 0.25	AMP 0.5	NEO 16	TYLT 32	TUL 32	TIL 16	CLI 0.5	ENRO 0.12
G	GEN 8	GEN 4	GEN 2	GEN 1	PEN 0.12	AMP 0.25	NEO 8	TYLT 16	TUL 16	TIL 8	CLI 0.25	POS
H	FFN 8	FFN 4	FFN 2	FFN 1	FFN 0.5	FFN 0.25	NEO 4	TYLT 8	TUL 8	TIL 4	POS	POS

ANTIMICROBICS

TIO	Ceftiofur
TIA	Tiamulin
CTET	Chlortetracycline
GEN	Gentamicin
FFN	Florfenicol
OXY	Oxytetracycline
PEN	Penicillin
AMP	Ampicillin
DANO	Danofloxacin
SDM	Sulphadimethoxine
NEO	Neomycin
SXT	Trimethoprim / sulfamethoxazo
SPE	Spectinomycin
TYLT	Tylosin tartrate
TUL	Tulathromycin
TIL	Tilmicosin
CLI	Clindamycin
ENRO	Enrofloxacin
POS	Positive Control