



SENSITITRE GRAM NEGATIVE NARMS PLATE FORMAT

Plate Code: **CMV1AGNF**

Plate Type: **MIC**

	1	2	3	4	5	6	7	8	9	10	11	12
A	FOX 32	AMI 32	CHL 16	AXO 64	AXO 0.25	CIP 2	GEN 16	NAL 16	TIO 2	FIS 32	KAN 64	AMP 2
B	FOX 16	AMI 16	CHL 8	AXO 32	AUG 32/16	CIP 1	GEN 8	NAL 8	TIO 1	FIS 16	KAN 32	AMP 1
C	FOX 8	AMI 8	CHL 4	AXO 16	AUG 16/8	CIP 0.5	GEN 4	NAL 4	TIO 0.5	SXT 4/76	KAN 16	STR 64
D	FOX 4	AMI 4	CHL 2	AXO 8	AUG 8/4	CIP 0.25	GEN 2	NAL 2	TIO 0.25	SXT 2/38	KAN 8	STR 32
E	FOX 2	AMI 2	TET 32	AXO 4	AUG 4/2	CIP 0.12	GEN 1	NAL 1	TIO 0.12	SXT 1/19	AMP 32	NEG
F	FOX 1	AMI 1	TET 16	AXO 2	AUG 2/1	CIP 0.06	GEN 0.5	NAL 0.5	FIS 256	SXT 0.5/9.5	AMP 16	POS
G	FOX 0.5	AMI 0.5	TET 8	AXO 1	AUG 1/0.5	CIP 0.03	GEN 0.25	TIO 8	FIS 128	SXT 0.25/4.75	AMP 8	POS
H	AMI 64	CHL 32	TET 4	AXO 0.5	CIP 4	CIP 0.015	NAL 32	TIO 4	FIS 64	SXT 0.12/2.38	AMP 4	POS

ANTIMICROBICS

AMI	Amikacin
AMP	Ampicillin
AUG	Amoxicillin/Clavulanic Acid
AXO	Ceftriaxone
CHL	Chloramphenicol
CIP	Ciprofloxacin
SXT	Trimethoprim/sulfamethoxazole
FOX	Cefoxitin
GEN	Gentamicin
KAN	Kanamycin
NAL	Nalidixic Acid
FIS	Sulfisoxazole
STR	Streptomycin
TET	Tetracycline
TIO	Ceftiofur
POS	Positive Control
NEG	Negative Control