

SENSITITRE CUSTOM PLATE FORMAT

Plate Code: **COMPAN1F**

	1	2	3	4	5	6	7	8	9	10	11	12
A	AMP 0.25	AMP 0.5	AMP 1	AMP 2	AMP 4	AMP 8	AMP 16	OXA+ 0.25	OXA+ 0.5	OXA+ 1	OXA+ 2	OXA+ 4
B	AUG2 4/2	AUG2 8/4	AUG2 16/8	AUG2 32/16	AMI 4	AMI 8	AMI 16	AMI 32	FOX 2	FOX 4	FOX 8	FOX 16
C	TIC 8	TIC 16	TIC 32	TIC 64	POD 2	POD 4	POD 8	POD 16	TIM2 8/2	TIM2 16/2	TIM2 32/2	TIM2 64/2
D	SXT 0.5/9.5	SXT 1/19	SXT 2/38	FOV 0.25	FOV 0.5	FOV 1	FOV 2	FOV 4	FOV 8	FAZ 4	FAZ 8	FAZ 16
E	GEN 1	GEN 2	GEN 4	GEN 8	IMI 1	IMI 2	IMI 4	IMI 8	CLI 0.5	CLI 1	CLI 2	CLI 4
F	PEN 0.06	PEN 0.12	PEN 0.25	PEN 0.5	PEN 1	PEN 2	PEN 4	PEN 8	DOX 2	DOX 4	DOX 8	POS
G	XNL 0.25	XNL 0.5	XNL 1	XNL 2	XNL 4	MAR 0.25	MAR 0.5	MAR 1	MAR 2	RIF 1	RIF 2	POS
H	ENRO 0.25	ENRO 0.5	ENRO 1	ENRO 2	ERY 0.5	ERY 1	ERY 2	ERY 4	CHL 4	CHL 8	CHL 16	POS

ANTIMICROBICS

AMP	Ampicillin
AUG2	Amoxicillin / clavulanic acid 2:1 ratio
TIC	Ticarcillin
SXT	Trimethoprim / sulfamethoxazole
GEN	Gentamicin
PEN	Penicillin
XNL	Ceftiofur
ENRO	Enrofloxacin
FOV	Cefovecin
AMI	Amikacin
POD	Cefpodoxime
IMI	Imipenem
ERY	Erythromycin
MAR	Marbofloxacin
OXA+	Oxacillin+2%NaCl
FOX	Cefoxitin
TIM2	Ticarcillin / clavulanic acid constant 2
CLI	Clindamycin
DOX	Doxycycline
CHL	Chloramphenicol
FAZ	Cefazolin
RIF	Rifampin
POS	Positive Control