

SENSITITRE CUSTOM PLATE FORMAT

2810101440

Plate Code: **COMPAN2F** Date: **24-Nov-10**

	1	2	3	4	5	6	7	8	9	10	11	12
A	AMP 0.12	AMP 0.25	AMP 0.5	AMP 1	CEP 2	CEP 4	CEP 8	OXA+ 0.25	OXA+ 0.5	OXA+ 1	OXA+ 2	OXA+ 4
B	AUG2 0.12/0.06	AUG2 0.25/0.12	AUG2 0.5/0.25	AUG2 1/0.5	AMI 4	AMI 8	AMI 16	AMI 32	FOX 2	FOX 4	FOX 8	FOX 16
C	TIC 8	TIC 16	TIC 32	TIC 64	POD 2	POD 4	POD 8	POD 16	TIM2 8/2	TIM2 16/2	TIM2 32/2	TIM2 64/2
D	SXT 0.5/9.5	SXT 1/19	SXT 2/38	FOV 0.25	FOV 0.5	FOV 1	FOV 2	FOV 4	FAZ 1	FAZ 2	FAZ 4	FAZ 8
E	GEN 1	GEN 2	GEN 4	GEN 8	IMI 1	IMI 2	IMI 4	IMI 8	CLI 0.5	CLI 1	CLI 2	CLI 4
F	PEN 0.06	PEN 0.12	PEN 0.25	PEN 0.5	PEN 1	PEN 2	PEN 4	PEN 8	DOX 2	DOX 4	DOX 8	POS
G	XNL 0.25	XNL 0.5	XNL 1	XNL 2	XNL 4	MAR 0.25	MAR 0.5	MAR 1	MAR 2	RIF 1	RIF 2	POS
H	ENRO 0.25	ENRO 0.5	ENRO 1	ENRO 2	ERY 0.5	ERY 1	ERY 2	ERY 4	CHL 4	CHL 8	CHL 16	POS

ANTIMICROBICS

AMP	Ampicillin
AUG2	Amoxicillin / clavulanic acid 2:1 ratio
TIC	Ticarcillin
SXT	Trimethoprim / sulfamethoxazole
GEN	Gentamicin
PEN	Penicillin
XNL	Ceftiofur
ENRO	Enrofloxacin
FOV	Cefovecin
CEP	Cephalothin
AMI	Amikacin
POD	Cefpodoxime
IMI	Imipenem
ERY	Erythromycin
MAR	Marbofloxacin
OXA+	Oxacillin+2%NaCl
FOX	Cefoxitin
TIM2	Ticarcillin / clavulanic acid constant 2
FAZ	Cefazolin
CLI	Clindamycin
DOX	Doxycycline
CHL	Chloramphenicol
RIF	Rifampin
POS	Positive Control

READ METHOD: ARIS/AUTOREAD/VIZION/SENSITOUCH/MANUAL

PRODUCT LABELLING - ORGANISM: Susceptibility testing of veterinary isolates

SUBSTRATE IN WELLS: GRAM NEGATIVE

RECONSTITUTION VOLUME: 50µl

INOCULUM CONCENTRATION: 1 x 10⁸ cfu/ml

PRODUCT LABELLING - USE: FOR VETERINARY USE ONLY

LABEL COLOUR: YELLOW

PLATE EXPIRY: 18 Months

PACK INSERT: To be assigned

ORGANISMS TO BE TESTED IN THE PLATE: GRAM POSITIVE AND GRAM NEGATIVE

MISCELLANEOUS: TREK INVENTORY

ALL SUSCEPTIBILITY CUSTOM PLATES HAVE A MINIMUM ORDER OF 500 PLATES. CUSTOM COMBO PLATES (ID + SUSCEPTIBILITY) HAVE A MINIMUM ORDER OF 2000 PLATES. A NON-CANCELLABLE PURCHASE ORDER MUST BE SUPPLIED AT THE TIME OF ORDERING. CUSTOM PLATES ORDERED WILL BE SHIPPED IN THEIR ENTIRETY TO THE ORDERING CUSTOMER IMMEDIATELY AFTER MANUFACTURE AND QC IS COMPLETED. PRODUCTION YIELD MAY VARY +/-10%. DELIVERY QUANTITY IS SUBJECT TO THE SAME VARIATION AND WILL BE INVOICED PRO-RATA. PAYMENT IN FULL IS DUE UPON RECEIPT OF THE PLATES.

CONTACT NAME: TELEPHONE NUMBER:

CUSTOMER SIGNATURE: DATE ORDERED