



SENSITITRE CAMPY PLATE FORMAT

Plate Code: **CAMPY**

Plate Type: **MIC**

	1	2	3	4	5	6	7	8	9	10	11	12
A	AZI 0.015	AZI 0.03	AZI 0.06	AZI 0.12	AZI 0.25	AZI 0.5	AZI 1	AZI 2	AZI 4	AZI 8	AZI 16	AZI 32
B	AZI 64	CIP 0.015	CIP 0.03	CIP 0.06	CIP 0.12	CIP 0.25	CIP 0.5	CIP 1	CIP 2	CIP 4	CIP 8	CIP 16
C	CIP 32	CIP 64	ERY 0.03	ERY 0.06	ERY 0.12	ERY 0.25	ERY 0.5	ERY 1	ERY 2	ERY 4	ERY 8	ERY 16
D	ERY 32	ERY 64	GEN 0.12	GEN 0.25	GEN 0.5	GEN 1	GEN 2	GEN 4	GEN 8	GEN 16	GEN 32	TET 0.06
E	TET 0.12	TET 0.25	TET 0.5	TET 1	TET 2	TET 4	TET 8	TET 16	TET 32	TET 64	FFN 0.03	FFN 0.06
F	FFN 0.12	FFN 0.25	FFN 0.5	FFN 1	FFN 2	FFN 4	FFN 8	FFN 16	FFN 32	FFN 64	NAL 4	NAL 8
G	NAL 16	NAL 32	NAL 64	TEL 0.015	TEL 0.03	TEL 0.06	TEL 0.12	TEL 0.25	TEL 0.5	TEL 1	TEL 2	TEL 4
H	TEL 8	CLI 0.03	CLI 0.06	CLI 0.12	CLI 0.25	CLI 0.5	CLI 1	CLI 2	CLI 4	CLI 8	CLI 16	POS

ANTIMICROBICS

AZI	Azithromycin
CIP	Ciprofloxacin
ERY	Erythromycin
GEN	Gentamicin
TET	Tetracycline
FFN	Florfenicol
NAL	Nalidixic Acid
TEL	Telithromycin
CLI	Clindamycin
POS	Positive Control

