

SENSITITRE CAMPYLOBACTER PLATE - EUCAMP

	1	2	3	4	5	6	7	8	9	10	11	12		
A	GEN 16	STR 1	STR 2	STR 4	STR 8	STR 16	GEN 16	STR 1	STR 2	STR 4	STR 8	STR 16		
B	GEN 8	CIP 4	TET 16	ERY 32	NAL 64	CHL 32	GEN 8	CIP 4	TET 16	ERY 32	NAL 64	CHL 32		
C	GEN 4	CIP 2	TET 8	ERY 16	NAL 32	CHL 16	GEN 4	CIP 2	TET 8	ERY 16	NAL 32	CHL 16		
D	GEN 2	CIP 1	TET 4	ERY 8	NAL 16	CHL 8	GEN 2	CIP 1	TET 4	ERY 8	NAL 16	CHL 8		
E	GEN 1	CIP 0.5	TET 2	ERY 4	NAL 8	CHL 4	GEN 1	CIP 0.5	TET 2	ERY 4	NAL 8	CHL 4		
F	GEN 0.5	CIP 0.25	TET 1	ERY 2	NAL 4	CHL 2	GEN 0.5	CIP 0.25	TET 1	ERY 2	NAL 4	CHL 2		
G	GEN 0.25	CIP 0.12	TET 0.5	ERY 1	NAL 2	POS CON	GEN 0.25	CIP 0.12	TET 0.5	ERY 1	NAL 2	POS CON		
H	GEN 0.12	CIP 0.06	TET 0.25	ERY 0.5	POS CON	POS CON	GEN 0.12	CIP 0.06	TET 0.25	ERY 0.5	POS CON	POS CON		

ANTIMICROBICS	
GEN	Gentamicin
STR	Streptomycin
CIP	Ciprofloxacin
TET	Tetracycline
ERY	Erythromycin
NAL	Nalidixic Acid
POS	Positive Control
CHL	Chloramphenicol