

SENSITITRE STAPHYLOCOCCI PLATE - EUST

	1	2	3	4	5	6	7	8	9	10	11	12
A	CLI 0.12	CLI 0.25	CLI 0.5	CLI 1	CLI 2	CLI 4	ERY 0.25	ERY 0.5	ERY 1	ERY 2	ERY 4	ERY 8
B	TET 0.5	TET 1	TET 2	TET 4	TET 8	TET 16	CIP 0.25	CIP 0.5	CIP 1	CIP 2	CIP 4	CIP 8
C	RIF 0.016	RIF 0.03	RIF 0.06	RIF 0.12	RIF 0.25	RIF 0.5	FOX 0.5	FOX 1	FOX 2	FOX 4	FOX 8	FOX 16
D	STR 4	STR 8	STR 16	STR 32	TIA 0.5	TIA 1	TIA 2	TIA 4	LZD 1	LZD 2	LZD 4	LZD 8
E	FUS 0.5	FUS 1	FUS 2	FUS 4	SYN 0.5	SYN 1	SYN 2	SYN 4	MUP 0.5	MUP 1	MUP 2	MUP 256
F	PEN 0.12	PEN 0.25	PEN 0.5	PEN 1	PEN 2	VAN 1	VAN 2	VAN 4	VAN 8	VAN 16	SMX 64	SMX 128
G	CHL 4	CHL 8	CHL 16	CHL 32	CHL 64	GEN 1	GEN 2	GEN 4	GEN 8	GEN 16	SMX 256	SMX 512
H	KAN 4	KAN 8	KAN 16	KAN 32	KAN 64	TMP 2	TMP 4	TMP 8	TMP 16	TMP 32	NEG CON	POS CON

ANTIMICROBICS

CLI	Clindamycin
TET	Tetracycline
RIF	Rifampin
STR	Streptomycin
FUS	Fusidate
PEN	Penicillin
CHL	Chloramphenicol
KAN	Kanamycin
TIA	Tiamulin
SYN	Quinupristin / dalfopristin
VAN	Vancomycin
GEN	Gentamicin
TMP	Trimethoprim
ERY	Erythromycin
CIP	Ciprofloxacin
FOX	Cefoxitin
LZD	Linezolid
MUP	Mupirocin
SMX	Sulfamethoxazole
NEG	Negative Control
POS	Positive Control