

SENSITITRE® CUSTOM PLATE FORMAT

Plate Code: **GN3F**

	1	2	3	4	5	6	7	8	9	10	11	12
A	AMI 64	A/S 32/16	FAZ 32	CEP 16	ETP 16	GEN 16	P/T4 128/4	SXT 4/76	TAZ 32	TGC 8	AXO 64	TET 16
B	AMI 32	A/S 16/8	FAZ 16	CEP 8	ETP 8	GEN 8	P/T4 64/4	SXT 2/38	TAZ 16	TGC 4	AXO 32	TET 8
C	AMI 16	A/S 8/4	FAZ 8	CEP 4	ETP 4	GEN 4	P/T4 32/4	SXT 1/19	TAZ 8	TGC 2	AXO 16	TET 4
D	AMI 8	A/S 4/2	FAZ 4	CEP 2	ETP 2	GEN 2	P/T4 16/4	SXT 0.5/9.5	TAZ 4	TGC 1	AXO 8	TET 2
E	AMP 32	AZT 32	FEP 32	MERO 8	FUR 32	CIP 4	FOX 32	POD 16	TAZ 2	TIM2 64/2	AXO 4	TET 1
F	AMP 16	AZT 16	FEP 16	MERO 4	FUR 16	CIP 2	FOX 16	POD 8	TAZ 1	TIM2 32/2	AXO 2	TET 0.5
G	AMP 8	AZT 8	FEP 8	MERO 2	FUR 8	CIP 1	FOX 8	POD 4	TOB 8	TIM2 16/2	AXO 1	NEG CON
H	AMP 4	AZT 4	FEP 4	MERO 1	FUR 4	CIP 0.5	FOX 4	POD 2	TOB 4	POS CON	POS CON	POS CON

ANTIMICROBICS

AMI	Amikacin
AMP	Ampicillin
A/S	Ampicillin/sulbactam
AZT	Aztreonam
FAZ	Cefazolin
FEP	Cefepime
CEP	Cephalothin
MERO	Meropenem
ETP	Ertapenem
FUR	Cefuroxime
GEN	Gentamicin
CIP	Ciprofloxacin
P/T4	Piperacillin / tazobactam constant 4
FOX	Cefoxitin
SXT	Trimethoprim / sulfamethoxazole
POD	Cefpodoxime
TAZ	Ceftazidime
TOB	Tobramycin
TGC	Tigecycline
TIM2	Ticarcillin / clavulanic acid constant 2
AXO	Ceftriaxone
TET	Tetracycline
NEG	Negative Control
POS	Positive Control

- Inoculum: 1 x 10⁵ cfu/ml
- Shelf Life: 24 mo.
- Ancillaries: Mueller-Hinton Broth w/ TES– 5 mL, 100/box (Part No. T3462), 10/box (Part No. T3462-10); Demineralized Water– 5ml, 100/box (Part No. T3339); 10/box (Part No. T3339-10)
- Plates per Box:10