

SENSITITRE® GRAM NEGATIVE ID PLATE FORMAT

Plate Code: **GNID**

Plate Type: **MIC**

	1	2	3	4	5	6	7	8	9	10	11	12
A	FR1	UREA	ORN	SORB	FR1	UREA	ORN	SORB	FR1	UREA	ORN	SORB
B	XYL	FR12	SUC	FR9	XYL	FR12	SUC	FR9	XYL	FR12	SUC	FR9
C	FR3	TREH	FR8	MANN	FR3	TREH	FR8	MANN	FR3	TREH	FR8	MANN
D	MALT	FR4	INOS	FR10	MALT	FR4	INOS	FR10	MALT	FR4	INOS	FR10
E	FR5	FRUC	ESCU	ATOL	FR5	FRUC	ESCU	ATOL	FR5	FRUC	ESCU	ATOL
F	ANOS	LYS	TDA	RAFF	ANOS	LYS	TDA	RAFF	ANOS	LYS	TDA	RAFF
G	FR7	ARG	FR6	CELL	FR7	ARG	FR6	CELL	FR7	ARG	FR6	CELL
H	MALN	PYR	CIT	AGMA	MALN	PYR	CIT	AGMA	MALN	PYR	CIT	AGMA

ACTIVE INGREDIENTS

FR1	Lysine AMC
UREA	Urea
ORN	Ornithine
SORB	Sorbitol
XYL	Xylose
FR12	2-Acetamido-2-deoxyglucoside MU + α -L-Arabinoside MU
SUC	Sucrose
FR9	β -D-glucoside MU
FR3	Phosphate MU
TREH	Trehalose
FR8	bis-Phosphate MU
MANN	Mannitol
MALT	Maltose
FR4	α -D-galactoside MU
INOS	Inositol
FR10	β -D-glucoside MU
FR5	Proline AMC
FRUC	Fructose
ESCU	Esculin
ATOL	Arabitol
ANOS	Arabinose
LYS	Lysine
TDA	TDA
RAFF	Raffinose
FR7	γ -Glutamine MU
ARG	Arginine
FR6	α -D-galactoside MU
CELL	Cellobiose
MALN	Malonate
PYR	Pyruvate
CIT	Citrate

- Inoculum: 1×10^5 cfu/ml
- Shelf Life: 12 mo.
- Ancillaries: Demineralized Water– 5ml, 100/box (Part No. T3339); 10/box (Part No. T3339-10); Mineral Oil– 250 mL (Part No. E10159)
- Plates per Box:10