

SENSITITRE® CUSTOM PLATE FORMAT

Plate Code: **GNXF**

	1	2	3	4	5	6	7	8	9	10	11	12
A	AMI 4	AZT 2	SXT 0.5/9.5	FEP 2	LEVO 1	CIP 0.25	MERO 1	DOR 0.12	DOR 0.25	DOR 0.5	DOR 1	DOR 2
B	AMI 8	AZT 4	SXT 1/19	FEP 4	LEVO 2	CIP .5	MERO 2	COL 0.25	COL 0.5	COL 1	COL 2	COL 4
C	AMI 16	AZT 8	SXT 2/38	FEP 8	LEVO 4	CIP 1	MERO 4	POL 0.25	POL 0.5	POL 1	POL 2	POL 4
D	AMI 32	AZT 16	SXT 4/76	FEP 16	LEVO 8	CIP 2	MERO 8	TAZ 1	TAZ 2	TAZ 4	TAZ 8	TAZ 16
E	TIM2 16/2	P/T4 8/4	GEN 1	TOB 2	MIN 2	DOX 2	FOT 1	FOT 2	FOT 4	FOT 8	FOT 16	FOT 32
F	TIM2 32/2	P/T4 16/4	GEN 2	TOB 4	MIN 4	DOX 4	TGC 0.25	TGC 0.5	TGC 1	TGC 2	TGC 4	TGC 8
G	TIM2 64/2	P/T4 32/4	GEN 4	TOB 8	MIN 8	DOX 8	ETP 1	ETP 2	ETP 4	ETP 8	ETP 16	POS
H	TIM2 128/2	P/T4 64/4	GEN 8	TOB 16	MIN 16	DOX 16	IMI 1	IMI 2	IMI 4	IMI 8	POS	POS

ANTIMICROBICS

AMI	Amikacin
TIM2	Ticarcillin / clavulanic acid constant 2
AZT	Aztreonam
P/T4	Piperacillin / tazobactam constant 4
SXT	Trimethoprim / sulfamethoxazole
MIN	Minocycline
FEP	Cefepime
GEN	Gentamicin
LEVO	Levofloxacin
TOB	Tobramycin
CIP	Ciprofloxacin
DOX	Doxycycline
MERO	Meropenem
FOT	Cefotaxime
TGC	Tigecycline
ETP	Ertapenem
IMI	Imipenem
DOR	Doripenem
COL	Colistin
POL	Polymixin B
TAZ	Ceftazidime
POS	Positive Control

- Inoculum: 1 x 10⁵ cfu/ml
- Shelf Life: 18 mo.
- Ancillaries: Mueller-Hinton Broth w/ TES– 5 ML, 100/box (Part No. T3462), 10/box (Part No. T3462-10); Demineralized Water– 5ml, 100/box (Part No. T3339); 10/box (Part No. T3339-10)
- Plates per Box: 10