

## SENSITITRE<sup>®</sup> HAEMOPHILLUS PLATE FORMAT

Plate Code: **HPB**

Plate Type: **MIC**

|          | 1            | 2           | 3            | 4           | 5           | 6           | 7                | 8            | 9           | 10           | 11          | 12          |
|----------|--------------|-------------|--------------|-------------|-------------|-------------|------------------|--------------|-------------|--------------|-------------|-------------|
| <b>A</b> | LEVO<br>4    | CLA<br>16   | PEN<br>1     | AXO<br>2    | AMP<br>4    | SPA<br>1    | SXT<br>2/38      | MERO<br>2    | FUR<br>8    | TET<br>4     | CHL<br>4    | AUG<br>16/8 |
| <b>B</b> | LEVO<br>2    | CLA<br>8    | PEN<br>0.5   | AXO<br>1    | AMP<br>2    | SPA<br>0.5  | SXT<br>1/19      | MERO<br>1    | FUR<br>4    | TET<br>2     | CHL<br>2    | AUG<br>8/4  |
| <b>C</b> | LEVO<br>1    | CLA<br>4    | PEN<br>0.25  | AXO<br>0.5  | AMP<br>1    | SPA<br>0.25 | SXT<br>0.5/9.5   | MERO<br>0.5  | FUR<br>2    | TET<br>1     | CHL<br>1    | AUG<br>4/2  |
| <b>D</b> | LEVO<br>0.5  | CLA<br>2    | PEN<br>0.12  | AXO<br>0.25 | AMP<br>0.5  | SPA<br>0.12 | SXT<br>0.25/4.75 | MERO<br>0.25 | FUR<br>1    | TET<br>0.5   | CHL<br>0.5  | AUG<br>2/1  |
| <b>E</b> | LEVO<br>0.25 | CLA<br>1    | PEN<br>0.06  | AXO<br>0.12 | AMP<br>0.25 | SPA<br>0.06 | SXT<br>0.12/2.38 | MERO<br>0.12 | FUR<br>0.5  | TET<br>0.25  | ERY<br>0.25 | ERY<br>0.5  |
| <b>F</b> | LEVO<br>0.12 | CLA<br>0.5  | PEN<br>0.03  | AXO<br>0.06 | AMP<br>0.12 | SPA<br>0.03 | SXT<br>0.06/1.19 | MERO<br>0.06 | FIX<br>0.12 | FIX<br>0.25  | FIX<br>0.5  | FIX<br>1    |
| <b>G</b> | LEVO<br>0.06 | CLA<br>0.25 | PEN<br>0.015 | AXO<br>0.03 | FEP<br>0.12 | FEP<br>0.25 | FEP<br>0.5       | FEP<br>1     | FEP<br>2    | A/S<br>1/0.5 | A/S<br>2/1  | NEG         |
| <b>H</b> | LEVO<br>0.03 | CLA<br>0.12 | FAC<br>4     | FAC<br>8    | FAC<br>16   | IMI<br>0.5  | IMI<br>1         | IMI<br>2     | IMI<br>4    | POS          | POS         | POS         |

### ANTIMICROBICS

|      |                               |
|------|-------------------------------|
| A/S  | Ampicillin/Sulbactam          |
| AMP  | Ampicillin                    |
| AUG  | Amoxicillin/Clavulanic Acid   |
| AXO  | Ceftriaxone                   |
| CHL  | Chloramphenicol               |
| CLA  | Clarithromycin                |
| SXT  | Trimethoprim/Sulfamethoxazole |
| ERY  | Erythromycin                  |
| FAC  | Cefaclor                      |
| FEP  | Cefepime                      |
| FIX  | Cefixime                      |
| FUR  | Cefuroxime                    |
| IMI  | Imipenem                      |
| LEVO | Levofloxacin                  |
| MERO | Meropenem                     |
| PEN  | Penicillin                    |
| SPA  | Sparfloxacin                  |
| TET  | Tetracycline                  |
| NEG  | Negative Control              |
| POS  | Positive Control              |

- Inoculum: 1 x 10<sup>5</sup> cfu/ml
- Shelf Life: 18 mo.
- Ancillaries: Mueller-Hinton Broth w/ TES- 5 mL, 100/box (Part No. T3462-05), 10/box (Part No. T3462-05-10); Mueller-Hinton Broth w/ LHB- 11 mL, 10/box (Part No. CP112-10); HTM Broth- 11 mL, 10/box (Part No. T3470)
- Plates per Box: 10