

SENSITITRE® HAEMOPHILUS PLATE FORMAT

Plate Code: **HPB**

Plate Type: **MIC**

	1	2	3	4	5	6	7	8	9	10	11	12
A	LEVO 4	CLA 16	PEN 1	AXO 2	AMP 4	SPA 1	SXT 2/38	MERO 2	FUR 8	TET 4	CHL 4	AUG 16/8
B	LEVO 2	CLA 8	PEN 0.5	AXO 1	AMP 2	SPA 0.5	SXT 1/19	MERO 1	FUR 4	TET 2	CHL 2	AUG 8/4
C	LEVO 1	CLA 4	PEN 0.25	AXO 0.5	AMP 1	SPA 0.25	SXT 0.5/9.5	MERO 0.5	FUR 2	TET 1	CHL 1	AUG 4/2
D	LEVO 0.5	CLA 2	PEN 0.12	AXO 0.25	AMP 0.5	SPA 0.12	SXT 0.25/4.75	MERO 0.25	FUR 1	TET 0.5	CHL 0.5	AUG 2/1
E	LEVO 0.25	CLA 1	PEN 0.06	AXO 0.12	AMP 0.25	SPA 0.06	SXT 0.12/2.38	MERO 0.12	FUR 0.5	TET 0.25	ERY 0.25	ERY 0.5
F	LEVO 0.12	CLA 0.5	PEN 0.03	AXO 0.06	AMP 0.12	SPA 0.03	SXT 0.06/1.19	MERO 0.06	FIX 0.12	FIX 0.25	FIX 0.5	FIX 1
G	LEVO 0.06	CLA 0.25	PEN 0.015	AXO 0.03	FEP 0.12	FEP 0.25	FEP 0.5	FEP 1	FEP 2	A/S 1/0.5	A/S 2/1	NEG
H	LEVO 0.03	CLA 0.12	FAC 4	FAC 8	FAC 16	IMI 0.5	IMI 1	IMI 2	IMI 4	POS	POS	POS

ANTIMICROBICS

A/S	Ampicillin/Sulbactam
AMP	Ampicillin
AUG	Amoxicillin/Clavulanic Acid
AXO	Ceftriaxone
CHL	Chloramphenicol
CLA	Clarithromycin
SXT	Trimethoprim/Sulfamethoxazole
ERY	Erythromycin
FAC	Cefaclor
FEP	Cefepime
FIX	Cefixime
FUR	Cefuroxime
IMI	Imipenem
LEVO	Levofloxacin
MERO	Meropenem
PEN	Penicillin
SPA	Sparfloxacin
TET	Tetracycline
NEG	Negative Control
POS	Positive Control

- Inoculum: 1 x 10⁵ cfu/ml
- Shelf Life: 18 mo.
- Ancillaries: Mueller-Hinton Broth w/ TES- 5 mL, 100/box (Part No. T3462-05), 10/box (Part No. T3462-05-10); Mueller-Hinton Broth w/ LHB- 11 mL, 10/box (Part No. CP112-10); HTM Broth- 11 mL, 10/box (Part No. T3470)
- Plates per Box: 10