



SENSITITRE MASTITIS PLATE FORMAT

Plate Code: **CMV1AMAF**

Plate Type: **MIC/BP**

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|----------|-------------|-------------|-------------|-------------|-----------|------------|-------------|-------------|-------------|-------------|-----------|------------|
| A | POS | POS | POS | PIRL 4 | TET 8 | XNL 4 | POS | POS | POS | PIRL 4 | TET 8 | XNL 4 |
| B | AMP 8 | PEN 8 | ERY 4 | PIRL 2 | TET 4 | XNL 2 | AMP 8 | PEN 8 | ERY 4 | PIRL 2 | TET 4 | XNL 2 |
| C | AMP 4 | PEN 4 | ERY 2 | PIRL 1 | TET 2 | XNL 1 | AMP 4 | PEN 4 | ERY 2 | PIRL 1 | TET 2 | XNL 1 |
| D | AMP 2 | PEN 2 | ERY 1 | PIRL 0.5 | TET 1 | XNL 0.5 | AMP 2 | PEN 2 | ERY 1 | PIRL 0.5 | TET 1 | XNL 0.5 |
| E | AMP 1 | PEN 1 | ERY 0.5 | P/N 8/16 | CEP 16 | SDM 256 | AMP 1 | PEN 1 | ERY 0.5 | P/N 8/16 | CEP 16 | SDM 256 |
| F | AMP 0.5 | PEN 0.5 | ERY 0.25 | P/N 4/8 | CEP 8 | SDM 128 | AMP 0.5 | PEN 0.5 | ERY 0.25 | P/N 4/8 | CEP 8 | SDM 128 |
| G | AMP 0.25 | PEN 0.25 | OXA+ 4 | P/N 2/4 | CEP 4 | SDM 64 | AMP 0.25 | PEN 0.25 | OXA+ 4 | P/N 2/4 | CEP 4 | SDM 64 |
| H | AMP 0.12 | PEN 0.12 | OXA+ 2 | P/N 1/2 | CEP 2 | SDM 32 | AMP 0.12 | PEN 0.12 | OXA+ 2 | P/N 1/2 | CEP 2 | SDM 32 |

ANTIMICROBICS

- AMP Ampicillin
- PEN Penicillin
- ERY Erythromycin
- OXA+ Oxacillin+2%NaCL
- PIRL Pirlimycin
- P/N Penicillin / novobiocin
- TET Tetracycline
- CEP Cephalothin
- XNL Ceftiofur
- SDM Sulphadimethoxine
- POS Positive Control

