

## SENSITITRE® CUSTOM PLATE FORMAT

Plate Code: **STP6F**

	1	2	3	4	5	6	7	8	9	10	11	12
<b>A</b>	MXF 1	MXF 2	MXF 4	MXF 8	PEN 0.03	PEN 0.06	PEN 0.12	PEN 0.25	PEN 0.5	PEN 1	PEN 2	PEN 4
<b>B</b>	LEVO 0.5	LEVO 1	LEVO 2	LEVO 4	MERO 0.25	MERO 0.5	MERO 1	MERO 2	AZI 0.25	AZI 0.5	AZI 1	AZI 2
<b>C</b>	TET 1	TET 2	TET 4	TET 8	ETP 0.5	ETP 1	ETP 2	ETP 4	ERY 0.25	ERY 0.5	ERY 1	ERY 2
<b>D</b>	FUR 0.5	FUR 1	FUR 2	FUR 4	AUG2 2/1	AUG2 4/2	AUG2 8/4	AUG2 16/8	SXT 0.5/9.5	SXT 1/19	SXT 2/38	SXT 4/76
<b>E</b>	AXO 0.12	AXO 0.25	AXO 0.5	AXO 1	AXO 2	LZD 0.25	LZD 0.5	LZD 1	LZD 2	LZD 4	VAN 0.5	VAN 1
<b>F</b>	FOT 0.12	FOT 0.25	FOT 0.5	FOT 1	FOT 2	FOT 4	CLI 0.12	CLI 0.25	CLI 0.5	CLI 1	VAN 2	VAN 4
<b>G</b>	DAP 0.06	DAP 0.12	DAP 0.25	DAP 0.5	DAP 1	DAP 2	FEP 0.5	FEP 1	FEP 2	FEP 4	FEP 8	POS
<b>H</b>	CHL 1	CHL 2	CHL 4	CHL 8	CHL 16	CHL 32	TGC 0.015	TGC 0.03	TGC 0.06	TGC 0.12	POS	POS

### ANTIMICROBICS

MXF	Moxifloxacin
LEVO	Levofloxacin
TET	Tetracycline
FUR	Cefuroxime
AXO	Ceftriaxone
FOT	Cefotaxime
DAP	Daptomycin
CHL	Chloramphenicol
PEN	Penicillin
MERO	Meropenem
ETP	Ertapenem
AUG2	Amoxicillin / clavulanic acid 2:1 ratio
LZD	Linezolid
CLI	Clindamycin
FEP	Cefepime
TGC	Tigecycline
AZI	Azithromycin
ERY	Erythromycin
SXT	Trimethoprim / sulfamethoxazole
VAN	Vancomycin
POS	Positive Control

- Inoculum: 1 x 10<sup>5</sup> cfu/ml
- Shelf Life: 18 mo.
- Ancillaries: Mueller-Hinton Broth w/ TES– 5 mL, 100/box (Part No. T3462), 10/box (Part No. T3462-10); Demineralized Water– 5ml, 100/box (Part No. T3339); 10/box (Part No. T3339-10)
- Plates per Box:10