



SENSITITRE URINARY PLATE FORMAT

Plate Code: **CMV1BURF**

Plate Code: **MIC**

	1	2	3	4	5	6	7	8	9	10	11	12
A	POS	POS	POS	AMP 256	AUG 256/128	ENRO 4	POS	POS	POS	AMP 256	AUG 256/128	ENRO 4
B	SXT 8/152	TET 128	LEX 256	AMP 128	AUG 128/64	ENRO 2	SXT 8/152	TET 128	LEX 256	AMP 128	AUG 128/64	ENRO 2
C	SXT 4/76	TET 64	LEX 128	AMP 64	AUG 64/32	ENRO 1	SXT 4/76	TET 64	LEX 128	AMP 64	AUG 64/32	ENRO 1
D	SXT 2/38	TET 32	LEX 64	AMP 32	AUG 32/16	ENRO 0.5	SXT 2/38	TET 32	LEX 64	AMP 32	AUG 32/16	ENRO 0.5
E	TIO 4	TET 16	LEX 32	AMP 16	AUG 16/8	ENRO 0.25	TIO 4	TET 16	LEX 32	AMP 16	AUG 16/8	ENRO 0.25
F	TIO 2	TET 8	LEX 16	AMP 8	AUG 8/4	ENRO 0.12	TIO 2	TET 8	LEX 16	AMP 8	AUG 8/4	ENRO 0.12
G	TIO 1	TET 4	LEX 8	AMP 4	AUG 4/2	ENRO 0.06	TIO 1	TET 4	LEX 8	AMP 4	AUG 4/2	ENRO 0.06
H	TIO 0.5	TET 2	LEX 4	AMP 2	AUG 2/1	ENRO 0.03	TIO 0.5	TET 2	LEX 4	AMP 2	AUG 2/1	ENRO 0.03

ANTIMICROBICS

- SXT Trimethoprim / sulfamethoxazole
- TIO Cefotiofur
- TET Tetracycline
- LEX Cephalixin
- AMP Ampicillin
- AUG Amoxicillin / clavulanic acid 2:1 ratio
- ENRO Enrofloxacin
- POS Positive Control