

# Magellan

B I O S C I E N C E S

October 7, 2010

Subject: New MIC plate (YO2IVD) replaces discontinued plate (YO2V)

Dear Customer,

We are proud to announce the availability of our newest Sensititre® susceptibility plate (Part No. YO2IVD) for antifungal susceptibility test procedures. This MIC plate is FDA cleared for visual read methods for yeast isolates. We have improved our plate by adding Caspofungin, making it more relevant to of our customers.

With the introduction of this new plate, we are discontinuing our previous MIC format, (Part No. YO2V) as we expect most of our users will replace it with the new format. Based on customer usage rates for Part No. YO2V, we estimate that we will be able to provide the previous format for an additional two month period to allow laboratories the opportunity to transition to the new format, complete necessary quality control procedures and manage additional internal changes. If for some reason, you do not wish to change, Sensititre plates can be customized to meet your specific needs. The new format, Part No. YO2IVD, will be made available to you at the same price as the YO2V product.

The YO2IVD two-isolate MIC plate contains the following antifungal agents and corresponding dilutions:

- Voriconazole: 0.008-4 µg/ml
- 5-Flucytosine 0.03- 32 µg/ml
- Fluconazole: 0.25-64 µg/ml
- Caspofungin: 0.015-8 µg/ml
- Itraconazole 0.03-2 µg/ml

We thank you for your continued use of Sensititre MIC plates.

If you have any additional questions regarding this change, please contact Technical Services at 1-800-642-7029. To process your standing order for the new or existing MIC format, please contact our Customer Service department in the U.S. at 1-800-871-8909, or internationally at +44 1342 318777.

Sincerely,

Jennifer Lorbach  
Marketing Director