

SENSITITRE[®] CUSTOM PLATE FORMAT

Plate Code: **YO2V**

	1	2	3	4	5	6	7	8	9	10	11	12
A	POS	VOR 0.008	VOR 0.015	VOR 0.03	VOR 0.06	VOR 0.12	VOR 0.25	VOR 0.5	VOR 1	VOR 2	VOR 4	VOR 8
B	FC 0.03	FC 0.06	FC 0.12	FC 0.25	FC 0.5	FC 1	FC 2	FC 4	FC 8	FC 16	FC 32	FC 64
C	FZ 0.12	FZ 0.25	FZ 0.5	FZ 1	FZ 2	FZ 4	FZ 8	FZ 16	FZ 32	FZ 64	FZ 128	FZ 256
D	IZ 0.008	IZ 0.016	IZ 0.03	IZ 0.06	IZ 0.12	IZ 0.25	IZ 0.5	IZ 1	IZ 2	IZ 4	IZ 8	IZ 16
E	POS	VOR 0.008	VOR 0.015	VOR 0.03	VOR 0.06	VOR 0.12	VOR 0.25	VOR 0.5	VOR 1	VOR 2	VOR 4	VOR 8
F	FC 0.03	FC 0.06	FC 0.12	FC 0.25	FC 0.5	FC 1	FC 2	FC 4	FC 8	FC 16	FC 32	FC 64
G	FZ 0.12	FZ 0.25	FZ 0.5	FZ 1	FZ 2	FZ 4	FZ 8	FZ 16	FZ 32	FZ 64	FZ 128	FZ 256
H	IZ 0.008	IZ 0.016	IZ 0.03	IZ 0.06	IZ 0.12	IZ 0.25	IZ 0.5	IZ 1	IZ 2	IZ 4	IZ 8	IZ 16

ANTIMICROBICS

POS	Positive Control
FC	5-Flucytosine
FZ	Fluconazole
IZ	Itraconazole
VOR	Voriconazole

- Inoculum: 1.5-8 x 10³ cfu/ml
- Shelf Life: 24 mo.
- Ancillaries: Demineralized Water– 5ml, 100/box (Part No. T3339), 10/box (Part No. T3339-10); YeastOne[®] Broth– 10/box (Part No. Y3462); YeastOne Interpretive Criteria Sheet (YO2-IT); YeastOne Quality Control Sheet (Part No. YO2-QC)
- Plates per Box:10