

SENSITITRE® CUSTOM PLATE FORMAT

Plate Code: **YO-9**

	1	2	3	4	5	6	7	8	9	10	11	12
A	POS	AND	AND	AND	AND	AND	AND	AND	AND	AND	AND	AB
		0.015	0.03	0.06	0.12	0.25	0.5	1	2	4	8	0.12
B	MF	MF	MF	MF	MF	MF	MF	MF	MF	MF	MF	AB
	0.008	0.015	0.03	0.06	0.12	0.25	0.5	1	2	4	8	0.25
C	CAS	CAS	CAS	CAS	CAS	CAS	CAS	CAS	CAS	CAS	CAS	AB
	0.008	0.015	0.03	0.06	0.12	0.25	0.5	1	2	4	8	0.5
D	FC	FC	FC	FC	FC	FC	FC	FC	FC	FC	FC	AB
	0.06	0.12	0.25	0.5	1	2	4	8	16	32	64	1
E	PZ	PZ	PZ	PZ	PZ	PZ	PZ	PZ	PZ	PZ	PZ	AB
	0.008	0.015	0.03	0.06	0.12	0.25	0.5	1	2	4	8	2
F	VOR	VOR	VOR	VOR	VOR	VOR	VOR	VOR	VOR	VOR	VOR	AB
	0.008	0.015	0.03	0.06	0.12	0.25	0.5	1	2	4	8	4
G	IZ	IZ	IZ	IZ	IZ	IZ	IZ	IZ	IZ	IZ	IZ	AB
	0.015	0.03	0.06	0.12	0.25	0.5	1	2	4	8	16	8
H	FZ	FZ	FZ	FZ	FZ	FZ	FZ	FZ	FZ	FZ	FZ	FZ
	0.12	0.25	0.5	1	2	4	8	16	32	64	128	256

ANTIMICROBICS

POS	Positive Control
AND	Anidulafungin
AB	Amphotericin B
MF	Micafungin
CAS	Caspofungin
FC	5-Flucytosine
PZ	Posaconazole
VOR	Voriconazole
IZ	Itraconazole
FZ	Fluconazole

- Inoculum: 1.5-8 x 10³ cfu/ml
- Shelf Life: 12 mo.
- Ancillaries: Demineralized Water– 5ml, 100/box (Part No. T3339); 10/box (Part No. T3339-10); YeastOne® Broth– 10/box (Part No. Y3462)
- Plates per Box:10