

VersaTREK[®]

Media Specifications

Blood Culture			
TREK Media	Catalog Number	Specimen Type	Specimen Volume
VersaTREK [™] REDOX1 [™] 80 ml with stir bar	7102-44	Blood Sterile body fluid	Range 0.1 to 10 ml
VersaTREK [™] REDOX2 [™] 80 ml	7103-44 No change	Blood Sterile body fluid	Range 0.1 to 10 ml
VersaTREK [™] REDOX 1 [™] EZ Draw [®] 40 ml with stir bar	7106-44	Blood Sterile body fluid	Range 0.1 to 5 ml
VersaTREK [™] REDOX 2 [™] EZ Draw [®] 40 ml	7107-44 No change	Blood Sterile body fluid	Range 0.1 to 5 ml
VersaTREK [™] REDOX 2M [™] EZ Draw [®] 40 ml	7110-44 No change	Blood Sterile body fluid	Range 0.1 to 5 ml



Mycobacteria			
TREK Media	Catalog Number	Specimen Type	Specimen Volume
VersaTREK [™] Myco	7111-42 No change	Blood Sterile body fluids Processed specimens	Up to 1 ml
VersaTREK [™] Myco GS	7112-42 No change	Myco Growth Supplement	N/A
VersaTREK [™] Myco PVNA	7113-42 No change	Myco Antibiotic Supplement	N/A
VersaTREK [™] Myco AS	7114-42 No change	Myco Antibiotic Supplement	N/A
VersaTREK [™] Myco Susceptibility Kit	7115-60 No change	Rifampin, Isoniazid, Ethambutol	Up to 1 ml
Myco PZA Kit*	7116-70 No change	Pyrazinamide	Up to 1 ml

*For research use only.
Not yet available for sale in the U.S.

VersaTREK® Automated Microbial Detection System



240 Model

Dimensions

52 in. width, 30 3/4 in. depth,
40 1/4 in. height
132 cm. width, 78 cm. depth,
103 cm. height

Unit Weight (Empty)

656 lbs.
304 kg.

Unit Weight (Fully Loaded)

766 lbs.
348 kg.

Capacity

96 to 240 bottles

528 Model

Dimensions

52 in. width, 30 3/4 in. depth,
76 1/4 in. height
132 cm. width, 78 cm. depth,
194 cm. height

Unit Weight (Empty)

1,212 lbs.
561 kg.

Unit Weight (Fully Loaded)

1,455 lbs.
660 kg.

Capacity

Up to 528 bottles



VersaTREK Windows® Software

VersaTREK Computer Package, with 15" Monitor 6110-31

TREK's Media Accessories

VersaTREK Myco Adaptors, 25/box	6040-30	Absorbent Sheet, each*	6444-30
Red Disable Location Cap, each*	6044-30	Label, Drawer, Numbered 1 - 22, each*	6452-30
VersaTREK Stylus, 3/pack*	6420-30	VersaTREK EZ-ID Labels, 100/pack	6450-35
Bottle Carry Rack, 24-position*	6424-30	Mounts, Thermometer, each*	6705-30
Thermometer, 18°C - 60°C, 125mm, each*	6440-30	VersaTREK Connectors, 50/box	7150-44

* These items are included in the VersaTREK Starter Kit (6499-31). Quantities of these items in the Starter Kit are: 5 red disable location caps, 1 VersaTREK stylus 3-pack, 1 bottle carry rack, 3 thermometers, 3 absorbent sheets, 1 set of drawer labels, and 6 thermometer mounts. One cable, CAT5E 25' Snagless Patch (6172-30), is also included with the Starter Kit. One complete Starter Kit is shipped with each VersaTREK instrument.

Specimen Collection, Preparation, and Safety

- VersaTREK REDOX 1 EZ Draw 40 ml with stir bar and REDOX 2 EZ Draw 40 ml provides the only FDA cleared, true, direct measured draw.
- VersaTREK Myco bottles are specially designed to resist breakage.
- Easy-to-operate with other collection devices.

Packaging and Storage Information

- All VersaTREK REDOX 80 ml and VersaTREK REDOX EZ Draw 40 ml bottles are packaged 50/box and should be stored at 15-30°C, away from light.
- VersaTREK Myco contains 50/box and should be stored at 15-30°C, away from light.
- VersaTREK Myco GS, Myco AS, and Myco PVNA are each packaged 5/box and should be stored 2-8°C. Myco GS must be stored in the dark.
- For technical information, refer to the technical inserts: VersaTREK REDOX, VersaTREK Myco, and VersaTREK Myco Susceptibility.



982 Keynote Circle, Suite 6
Cleveland, Ohio 44131
800.871.8909
FAX: 216.351.5456

Imberhorne Lane, East Grinstead
West Sussex, RH19 1QX England
44-1342-318777
FAX: 44-1342-318666

www.trekds.com